

2025
ESSENTIAL
INSTRUMENTATION
CATALOG




TABLE OF CONTENTS



SECTION A DIAGNOSTIC - 4



SECTION B PERIODONTAL - 16



SECTION C PERIODONTAL SURGERY - 32




SECTION D IMPLANT & SURGICAL - 44




SECTION E ELEVATORS & LUXATORS - 112



SECTION F FORCEPS - 132



SECTION G OPERATIVE - 150



SECTION H ENDODONTIC - 176




SECTION I CROWN INSTRUMENTS - 186



SECTION J ORTHODONTIC - 196



SECTION K GENERAL INSTRUMENTS - 214



SECTION L KITS & SETS - 228

HANDLE OPTIONS

SINGLE-END



1EZ S/E (1/4" 6mm)
STAINLESS STEEL



4EZ S/E (1/4" 6mm)
STAINLESS STEEL



6EZ S/E (1/4" 6mm)
STAINLESS STEEL

DOUBLE-END



3EZ D/E (1/4" 6mm)
STAINLESS STEEL



4EZ D/E (5/16" 8mm)
STAINLESS STEEL



6EZ D/E (3/8" 10mm)
STAINLESS STEEL

ELEVATOR



SMALL ELEVATOR
STAINLESS STEEL



LARGE ELEVATOR
STAINLESS STEEL

We use the highest quality stainless steel for our handles. Our handles are guaranteed against tarnishing, cracking and deforming. **Handle options are subject to change due to availability.**



DIAGNOSTIC



A

EXPLORERS

Explorers are used to detect & explore caries, abrasions, calculus, pocket characteristics, furcation, and restoration. Moderate tips are used for disease detection and hard tips are used for removing calculus & excessive cement.



1005 EXP23 1EZ



1005E EXP23 6EZ



1009
EXP3A
3EZ

1009E
EXP3A
6EZ

1010
EXP6
3EZ

1010E
EXP6
6EZ

1011
EXP6A
3EZ

1011E
EXP6A
6EZ

1012
EXP54
3EZ

1012E
EXP54
6EZ

DIAGNOSTIC | EXPLORERS



1006 EXD5 3EZ



1006E EXD5 6EZ



1006FS EXP X POWER FLEX EXD 5 FS (5/PKG)




1006L EXD5L 3EZ THICKER VERSION



1015 EXD2 3EZ



1007 EXP TU17/23 3EZ



1015E EXD2 6EZ



1007E EXP TU17/23 6EZ



1017 EXD3 3EZ




1013 EXD1



1017E EXD3 6EZ



1014 EXD2A 3EZ



1018 ODU 11/12 3EZ



1014E EXD2A 6EZ



1018E ODU 11/12 6EZ

DIAGNOSTIC | EXPLORERS




1019 EXD 11/12AF 3EZ



1022E EXD 6 6EZ



1019E EXD 11/12AF 6EZ



1023 EXD 23/6 3EZ



1020 EXD3ES 3EZ



1023E EXD 23/6 6EZ



1020E EXD3ES 6EZ




1026 DG 16 3EZ



1021 EXD 3CH 3EZ



1026E DG 16 6EZ



1021E EXD 3CH 6EZ



1029 23/DG16 3EZ



1022 EXD 6 3EZ

COLOR PROBES

Probes are used to check the depth of the periodontal pocket either vertically or horizontally. A length of 6mm is the standard control point for assessing a patient's periodontal condition. A length of less than 6mm indicates that the periodontal condition is mild; a length of more than 6mm indicates that it is critical.



1101 CP8 (3-6-8-11) 3EZ



1101E CP8 (3-6-8-11) 6EZ



1102
CP10 3EZ
(1-2-3-5-7-8-9-10mm)



1102E
CP10 6EZ
(1-2-3-5-7-8-9-10mm)



1103
CP11 3EZ
(3-6-8-11mm)



1103E
CP11 6EZ
(3-6-8-11mm)



1104
CP12 3EZ
(3-6-9-12mm)



1104B
CP12 X
(3-6-9-12mm)



1104PS
CP12 5PK
(3-6-9-12mm)

DIAGNOSTIC | COLOR PROBES



1104E
CP12 6EZ
(3-6-9-12mm)



1105
CP15UNC 3EZ
(1-2-3-4-5-6-7-8-9-10-11-12-13-14-15mm)



1105BN
CP15UNC
(1-2-3-4-5-6-7-8-9-10-11-12-13-14-15mm)



1105E
CP15UNC 6EZ
(1-2-3-4-5-6-7-8-9-10-11-12-13-14-15mm)



1106
CP0 3EZ
(3-6-8mm)



1106E
CP0 6EZ
(3-6-8mm)



1107
CPQW 3EZ
(1-2-3-5-7-8-9-10mm)



1107E
CPQW 6EZ
(1-2-3-5-7-8-9-10mm)



1108
CP11.5B 3EZ
(3.5-5.5-8.5-11.5mm)



1108E
CP11.5B 6EZ
(3.5-5.5-8.5-11.5mm)



1109
CPN11.5 3EZ
(3.5-5.5-8.5-11.5mm)



1109E
CPN11.5 3EZ
(3.5-5.5-8.5-11.5mm)



1113
CPGF 3EZ
(1-2-3-5-7-8-9-10mm)



1113E
CPGF 6EZ
(1-2-3-5-7-8-9-10mm)



1137
CP2 3EZ
(1-2-4-6-8-10-12mm)



1137E
CP2 6EZ
(1-2-4-6-8-10-12mm)

EXPLORERS/ PROBES

Explorers/Probes are used to check the depth of the periodontal pocket either vertically or horizontally. A length of 6mm is the standard control point for assessing a patient's periodontal condition. A length of less than 6mm indicates that the periodontal condition is mild; a length of more than 6mm indicates that it is critical.



1133E EXP23/UNC15 6EZ
(1-2-4-5-6-7-8-9-10-11-12-13-14mm)



Goldman-Fox

1120E GF/CPQW 6EZ
(1-2-3-5-7-8-9-10mm)

Williams




1127E EXP23/CPQW 6EZ
(1-2-3-5-7-8-9-10mm)




1124 PR2N, NON-COLOR 3EZ



1128 EXP23/CP8 3EZ
(3-6-8-11mm)



1124E PR2N, NON-COLOR 6EZ



1128E EXP23/CP8 6EZ
(3-6-8-11mm)



1125 CPQN 3EZ
(3-6-9-12mm)



1129 EXP23/CP10 3EZ
(1-2-3-5-7-8-9-10mm)




1125E CPQN 6EZ
(3-6-9-12mm)



1129E EXP23/CP10 6EZ
(1-2-3-5-7-8-9-10mm)

DIAGNOSTIC | EXPLORERS/PROBES



1151 X PROBE SET OF 3



MIRRORS/ MIRROR HANDLES

Enhance clarity with our Power Mirrors. Coming in a variety of sizes, styles, and handles, improve efficiency through your view.



6908A OCTAGON C/S WITH RULER



6908E ENDO C/S WITH MEASUREMENT



6910 4EZ C/S



6911 6EZ C/S



6911BE X IMPLANT SURGERY



6912 6EZ D/E DOUBLE-SIDED



6913 X POWER C/S SURGICAL



6933 LARGE ELEVATOR HANDLE



6950B MIRROR #4 22.8mm C/S
FRONT SURFACE BOX (12 EA)



6951B MIRROR #5 24.8mm C/S
FRONT SURFACE BOX (12 EA)



6951DB MIRROR #5 24.8mm C/S
DOUBLE-SIDED BOX (6 EA)



6953 MIRROR X POWER SPECIAL
SURGICAL LONG SHAPE C/S
DOUBLE-SIDED 22x37mm (6EA)

DIAGNOSTIC | MIRRORS/MIRROR HANDLES



6960 #5 MIRROR WITH S/E THIN RULER MIRROR HANDLE (6908A & 6951)



6961 MIRROR #5 WITH 4EZ HANDLE (6951+6910)



6962 MIRROR #5 WITH 6EZ HANDLE (6951+6911)




6963 #4 MIRROR C/S FRONT SURFACE & #5 MIRROR C/S DOUBLE SURFACE WITH D/E MIRROR HANDLE (6911D, 6951, 6951D)



6964 #5 MIRROR DOUBLE-SIDED WITH 6EZ S/E MIRROR HANDLE (6911 & 6951D)



6969 #4 MIRROR SINGLE-SIDED WITH S/E THIN RULER MIRROR HANDLE (6908A & 6950)



6971 POWER MIRROR WITH HANDLE (6953 & 6911BE)



6931 PARA MIRROR ANTERIOR 5mm MEASUREMENT LINE STAINLESS



6932 PARA MIRROR POSTERIOR 5mm MEASUREMENT LINE STAINLESS



PERIODONTAL



B

SICKLE & JACQUETTE SCALERS

Primarily used in the Anterior region, get the deep cleaning and bite you've been craving with our Power Sickle & Jacquette Scalers.

ANTERIOR



1303 N128



1201 H6/7




1201B H6/7
BLACK EDITION HANDLE SS TIPS



1201BT H6/7
BLACK EDITION HANDLE TITANIUM TIPS



1201T H6/7
SS HANDLE TITANIUM TIPS




1202 H5/U15



1203 H5/33



1203B H5/33
BLACK EDITION HANDLE SS TIPS




1203T H5/33
SS HANDLE TITANIUM TIPS




1205 U15/30



1206 H5/33




1213 J30/33



1223 N137



1223B N137 BLACK EDITION



1224 MO/00





1225 NEVI-1/H5




1225M NEVI-1/H5 MINI

PERIODONTAL | SICKLE & JACQUETTE SCALERS



POSTERIOR




SCALERS & CURETTES

Our scalers are built with a pointed tip to grab onto supra-gingival surfaces and in tight interproximal spaces. Our curettes have a rounded toe which allow for scaling within a periodontal pocket (sub-gingivally) or supra-gingivally depending on the design of curette.


BALL SCALERS




1313 B2/3




1314 CK6



1316 B5/6



1315 C12/3



1317 CK6D

McCALLS



1322 MC13/14 ROUNDED (AVAILABLE POINTED 1322S)



1323 MC17/18 ROUNDED (AVAILABLE POINTED 1326PT)



1322ST MC13/14S TITANIUM POINTED



1325 MTC13/14 (AVAILABLE POINTED 1325S)

SINGLE-ENDED SCALERS



GOLDMAN FOX



CATTONI



UNIVERSAL CURETTES

Remove calculus and access everywhere in the mouth with our precision universal cures. Versatility in design allows for access across all 32 teeth.



1501 IU13/14




1502 IU17/18



1512 ME1



1503 BUNT1/2



1506 BT13/14S




1504 BUNT5/6



1507 RULE3/4



1505 YG7/8




1508 411/412

PERIODONTAL | UNIVERSAL CURETTES


LANGER




1514 LANGER1/2



1515 LANGER3/4




1516 LANGER5/6




1517 LANGER17/18

LANGER AFTER FIVES



1522 LANGER1/2 AF



1523 LANGER3/4 AF




1524 LANGER5/6 AF

COLUMBIA



1414 COL2R/2L



1415 COL4R/4L



1415T COL4R/4L TITANIUM TIPS



1416 COL13/14



1416C COL13/14C




1416T COL13/14T TITANIUM TIPS

BARNHART



1421 BH1/2



1422 BH5/6



1422T BH5/6 TITANIUM TIPS

GRACEY CURETTES

Gracey curettes are designed to improve access to the distal surfaces of posterior teeth and used for the removal of both supragingival and subgingival calculus in specific areas of the mouth including deep scaling, root planing, and periodontal debridement.


PROPHY




1518 PROPHE GR1/2



1519 PROPHE GR3/4



1520 PROPHE GR5/6



1521 PROPHE GR7/8




1601
GR1/2


1610
AF1/2

1601S
SAS1/2


GRACEY



1601 GR1/2




1607 GR13/14




1602 GR3/4




1608 GR15/16



1603 GR5/6



1609 GR17/18



1604 GR7/8



1623 GR11/14



1605 GR9/10




1624 GR12/13



1606 GR11/12

GRACEY MINI



1601M GR1/2 MINI



1607M GR13/14 MINI



1602M GR3/4 MINI



1608M GR15/16 MINI



1603M GR5M/6M MINI



1609M GR17/18 MINI



1604M GR7/8 MINI



1623M GR11/14 MINI



1605M GR9/10 MINI



1624M GR12/13 MINI




1606M GR11M/12M MINI

GRACEY RIGID



1601R GR1/2 RIGID



1607R GR13/14 RIGID



1602R GR3/4 RIGID



1608R GR15/16 RIGID



1603R GR5/6 RIGID



1609R GR17/18 RIGID



1604R GR7/8 RIGID



1623R GR11/14 RIGID



1605R GR9/10 RIGID




1624R GR12/13 RIGID



1606R GR11/12 RIGID

AFTER FIVE GRACEY



1610 AFG1/2



1615 AFG13/14



1611 AFG3/4




1616 AFG15/16



1612 AFG5/6




1617 AFG11/14



1613 AFG7/8



1618 AFG12/13




1614 AFG11/12

AFTER FIVE GRACEY RIGID




1610R AFG1/2 RIGID



1615R AFG13/14 RIGID



1611R AFG3/4 RIGID



1616R AFG15/16 RIGID



1612R AFG5/6 RIGID




1613R AFG7/8 RIGID



1617R AFG11/14 RIGID




1614R AFG11/12 RIGID




1618R AFG12/13 RIGID

SAS GRACEY




1610S SAS1/2




1615S SAS13/14



1611S SAS3/4



1616S SAS15/16



1612S SAS5/6



1613S SAS7/8



1617S SAS11/14



1614S SAS11/12



1618S SAS12/13

SASR GRACEY RIGID



1610SR SAS1/2 RIGID



1615SR SAS13/14 RIGID



1611SR SAS3/4 RIGID



1616SR SAS15/16 RIGID



1612SR SAS5/6 RIGID




1617SR SAS11/14 RIGID



1613SR SAS7/8 RIGID




1614SR SAS11/12 RIGID




1618SR SAS12/13 RIGID


SUB-O GRACEY




1619 SUB-O



1621 SUB11/12



1620 SUB1/2




1622 SUB13/14

MICRO MINI GRACEY



1625 SMS 1/2




1627 SMS 11/12



1626 SMS 7/8



1628 SMS 13/14



DEBRIDEMENT & FURCATION


O'Heir cures are used for removing residual calculus after ultrasonic scaling for deep pockets. Furcation cures are used for curetting the infected areas of periodontal disease.



1701 O'HEHIR DC1/2




1702 O'HEHIR DC3/4




1703 O'HEHIR DC5/6



1704 O'HEHIR DC7/8



1705 OR6/7



1706 H8/9



1707 U5/TG




1708 U5/13K



1711 DM1




1712 DM2




1721 DP5/6




1722 DP9/11




1723 HF3/7




1724 OR10/11



1725 OR12/13



1726 HF5/11



1727 HF5/11



1709
XS15



1710
7K



1718
CU62



1720
ZER

SHARPENING STONES

The coarse grained India stone is for the initial grinding of very blunt instruments. The shape of the stone enables the grinding of most of the surface and radii of the instrument. Fine grinding should then be carried out with an Arkansas stone.



7423

SHARPENING STONE INDIA SLIP



7424

SHARPENING STONE ARKANSAS FLAT



PERIODONTAL SURGERY



C

SURGICAL & KNIVES

Periodontal dental surgery knives are used for a variety of procedures, including creating bevel incisions, contouring the gingiva, excising interproximal tissue, removing excess gum tissue, trimming gum shape, and promoting periodontal tissue regeneration and repair.



1801 K15/16



1802 KGF8




1806 BUCK5/6



1803 KGF11




1807 KGF7



1804 KSA1/2



1808 KGF9



1805 KUSC19/20



1809 BUCK3/4



1810 KPO1/2



1811 SOLT1/2



1812 KOR1/2



1812L KOR1/2 LONG BLADE



1812M KOR1/2 MICRO BLADE



1813 KCCKE



1816D KM3/4

PERIODONTAL FILES

Periodontal files are dental instruments used to break up calculus deposits and other dental work including debridement in deep pockets and scaling and root planing. They can also be used to remove amalgam overhangs and overhanging restorations.



1818 F9/10S



1819 F11/12S



1820 F1/2S



1821 F3/4S

SURGICAL CURETTES

Surgical curettes are used to remove subgingival calculus, plaque, and diseased tissue from the teeth and periodontal pockets during periodontal surgery procedures. They can also be used to remove unwanted debris or cysts from the tooth socket, transfer bone grafts, and clean out infected cavity areas.



1901 PR1/2



1902 KRA1



1903 KRA2



1904 KRA3



1905 KN3/4



1907 KIR4/5



1908 KIR8/9



1918 SS3/6



1919 SS9/10



1920 SS11/12

PERIOTOMES

A perioste is a dental instrument used to remove a tooth by cutting the periodontal ligament (PDL) that holds it in place. The serrated blade cuts the PDL easily, resulting in minimal tissue damage.




1923 PT2



1923T PT2 TITANIUM



1924T PT1 TITANIUM



1924 PT1



1924TF PT1 GOLD TITANIUM



1925
PT3



1925T
PT3
TITANIUM



1926T
PT4
TITANIUM



1927
PT5

PERIODONTAL SURGERY | PERIOTOMES



1927T PT5 TITANIUM



1928 MS1



1929A ELIO-1



1929B ELIO-2



2043 PPael



2043T PPael TITANIUM



2044 PPaelA



2044T PPaelA TITANIUM

PERIODONTAL SURGERY | PERIOTOMES

1930FSG

X PERIOTOME POWER FLEX & ROOT EXTRACTION SET OF 8



1930-1G POWER FLEX PERIOTOME STR S/E



1930-2G POWER FLEX PERIOTOME CVD S/E



1930-3G X ROOT EXTRACTION TITANIUM SCREW ANTERIOR



1930-4G POWER FLEX PERIOTOME FPT-2 D/E



1930-5G POWER FLEX PERIOTOME FPT-1 D/E



1930-6G POWER FLEX PERIOTOME TWIST D/E



1930-3SG X ROOT EXTRACTION TITANIUM SCREW POSTERIOR




7821G X ROOT EXTRACTION TITANIUM SCREW HANDPIECE ANCHOR 30mm

4032

POWER TWIST COMPLETE SET OF 6




4024 POWER TWIST ELEVATOR STR-R



4023 POWER TWIST ELEVATOR STR-L



3897G POWER TWIST PERIOTOME 4mm POSTERIOR, CVD-R



3896G POWER TWIST PERIOTOME 4mm POSTERIOR, CVD-L



3895G POWER TWIST PERIOTOME 4mm ANTERIOR, STR-R




3894G POWER TWIST PERIOTOME 4mm ANTERIOR, STR-L

SURGICAL CHISELS

Surgical chisels are used to split, remove and reshape bone, reflect flaps, and remove secondary palatal flaps. They can also be used to elevate fractured root tips and remove remaining fibers on teeth as well as slightly remodeling bone margins, removing soft tissue, and for significant deposits of supragingival calculus.




2001 13K/TG 3.5/3mm




2001T 13K/TG 3.5/3mm TITANIUM



2002T 3K/LR 4.5/3.5mm TITANIUM



2002 13K/LR 4.5/3.5mm




2003 KN1/2




2002E 13K/LR 4.5/3.5mm BLACK EDITION



2004 C36/37



2002L 13K/13KL 4/6mm



2005 44/45

PERIODONTAL SURGERY | SURGICAL CHISELS




2006 TG




2007 SOLT3/4



2008 TG-O



2010 OCH1



2011 OCH2



2014 OCH4



2015 FEDI1



2016 FEDI2



IMPLANT & SURGICAL



D

BLACK TITANIUM INSTRUMENTS

Our black titanium instruments are enhanced in strength and durability up to three times more than ordinary tips. They are easy to clean and require less sharpening than other instruments. The titanium coating provides a non-stick surface.



1601T GR1/2



1201BT H6/7 BLACK EDITION HANDLE




1209T 204SD




1201T H6/7



1226T NEVI-2



1203T H5/33




1322ST MC13/14S POINTED



1208BT 204S BLACK EDITION HANDLE



1415BT COL4R/4L BLACK HANDLE



1208T 204S




1415T COL4R/4L


IMPLANT & SURGICAL | BLACK TITANIUM INSTRUMENTS



1416BT COL13/14 BLACK HANDLE



1608T GR15/16




1416T COL13/14T



1609T GR17/18




1422T BH5/6



1610T AFG1/2



1601BT GR1/2 BLACK HANDLE



1614T AFG11/12




1602T GR3/4




1615T AFG13/14




1606BT GR11/12 BLACK HANDLE



1620T SUB1/2




1606T GR11/12




1621T SUB11/12



1607BT GR13/14 BLACK HANDLE



1607T GR13/14




1622T SUB13/14

SURGICAL CHISELS


Surgical chisels are used to prepare the bone for dental implants by contouring, smoothing, or removing the bone so that implant material can easily be inserted. They can also be used to split teeth, flatten enamel, and remove secondary palatal flaps.




2001 13K/TG 3.5/3mm




2001T 13K/TG 3.5/3mm TITANIUM



2002 13K/LR 4.5/3.5mm



2002L 13K/13KL 4/6mm




2002T 3K/LR 4.5/3.5mm TITANIUM



2002E 13K/LR 4.5/3.5mm BLACK EDITION



2010 OCH1



2011 OCH2

BLOCK FORCEPS



2067
FORCEP BONE BLOCK
19cm STR




2067N
FORCEP BONE BLOCK
19cm CVD




2068
FORCEP CORTICAL BONE
BLOCK 19cm CVD

OSTEOTOMES/ BONE SHAPING


An osteotome is a surgical instrument used to cut, compress, or deform bone during dental implant procedures or bone grafting. They can also be used to create implant sites and elevate the sinus floor.




2061 OSTEOTOME 6.0mm JOVANOVIC




2064 RIDGE SPLIT CHISEL 3mm SINGLE/CUT (8-10-13-16mm)




2064D RIDGE SPLIT CHISEL 3mm DOUBLE/CUT (8-10-13-16mm)




2065 RIDGE SPLIT CHISEL 6mm SINGLE/CUT (8-10-13-16mm)



2065D RIDGE SPLIT CHISEL 6mm DOUBLE/CUT (8-10-13-16mm)



2066 RIDGE SPLIT CHISEL 9mm SINGLE/CUT (8-10-13-16mm)




2066D RIDGE SPLIT CHISEL 9mm DOUBLE/CUT (8-10-13-16mm)

OSTEOTOMES/ SINUS LIFT


An osteotome is a surgical instrument used to cut, compress, or deform bone during dental implant procedures or bone grafting. They can also be used to create implant sites and elevate the sinus floor.




2057 OSTEOTOME 4mm STR (7-10-13-15-18mm)



2058 OSTEOTOME 4mm CVD (7-10-13-15-18mm)



2059 OSTEOTOME 7.5mm STR (7-10-13-15-18mm)



2060 OSTEOTOME 7.5mm CVD (7-10-13-15-18mm)

IMPLANT & SURGICAL | OSTEOTOMES



2050 OSTEOTOME 3.2mm CVD (8-10-13-15-18mm)



2051 OSTEOTOME 4.0mm CVD (8-10-13-15-18mm)



2052 OSTEOTOME 1.8mm CVD (8-10-13-15-18mm)




2053 OSTEOTOME 2.5mm CVD (8-10-13-15-18mm)



2054 OSTEOTOME 3.0mm CVD (8-10-13-15-18mm)



2055 OSTEOTOME 4.0mm CVD (8-10-13-15-18mm)



2056 OSTEOTOME R0.38mm CVD (8-10-13-15-18mm)

IMPLANT & SURGICAL | OSTEOTOMES



2050S OSTEOTOME 3.2mm STR (8-10-13-15-18mm)



2051S OSTEOTOME 4.0mm STR (8-10-13-15-18mm)



2052S OSTEOTOME 1.8mm STR (8-10-13-15-18mm)



2053S OSTEOTOME 2.5mm STR (8-10-13-15-18mm)



2054S OSTEOTOME 3.0mm STR (8-10-13-15-18mm)



2055S OSTEOTOME 4.0mm STR (8-10-13-15-18mm)



2056S OSTEOTOME R0.38mm STR (8-10-13-15-18mm)

SINUS LIFT INSTRUMENTS

Sinus lift instruments are utilized for gently lifting the sinus membrane away from the floor and lateral portion of the maxillary sinus. Coming in a variety of shapes and sizes, this ensures clean lift at any angle necessary.



2021 M6577



2023 SINUS LIFT 1



2024 SINUS LIFT 2



2025 SINUS LIFT 3



2026 SINUS LIFT 4

IMPLANT & SURGICAL | SINUS LIFT INSTRUMENTS



2027 SINUS LIFT 5



2062 SINUS LIFT 6




2063 SINUS LIFT 7



2069 MUSHROOM 2.8mm/3.8mm (2-4-5-6-8-10)



2070 MUSHROOM 4.8mm/5.8mm (2-4-5-6-8-10)



2037T BONE GAUGE 4-6-8-10-12-14-16mm TITANIUM

2070S

SINUS LIFT MEMBRANE ELEVATION KIT SET OF 7 PCS



1970 SINUS BONE GRAFT 3.87SP/2.8mm (2-4-5-6-8-10mm)

2062 SINUS LIFT 6

2063 SINUS LIFT 7

2069 MUSHROOM 2.8mm/3.8mm (2-4-5-6-8-10)

2035T BONE CARRIER TITANIUM

2038T BONE GRAFT/CARRIER 2-4-6-8-10mm TITANIUM

BONE GRAFT PACKERS

Dental bone graft packers are specialized tools used in dental surgeries to assist with placing and compressing bone graft material into a surgical site.



2032 BONE PLUGGER IMPLANT LABIANCA 3X1.5mm



2033 BONE PLUGGER IMPLANT JOVANOVIC 3mm



2033A BONE GRAFT PACKER 1.0/1.4mm (8-11-13-15-18mm) TITANIUM



2033B BONE GRAFT PACKER 1.7/2.0mm (8-11-13-15-18mm) TITANIUM



2033C BONE GRAFT PACKER 2.0/2.7mm (8-11-13-15-18mm) TITANIUM

IMPLANT & SURGICAL | BONE GRAFT PACKERS



2033N BONE GRAFT PACKER 3.3/4mm (8-11-13-15-18mm)



2033NT BONE GRAFT PACKER 3.3/4mm TITANIUM



2033NTS BONE GRAFT PACKER 3.3/4mm TITANIUM SERRATED



2034S BONE GRAFT PACKER 3/5mm SERRATED



2035T BONE CARRIER TITANIUM



2036S BONE GRAFT PACKER 4/6mm SERRATED



2038T BONE GRAFT/CARRIER 2-4-6-8-10mm TITANIUM

BONE INSTRUMENTS

Bone instruments are essential tools used in oral and maxillofacial surgeries, implantology, and periodontal treatments to manage, modify, and augment bone.




2071 KERRISON FORCEP 3mm
2072 KERRISON FORCEP 4mm



2030 BONE SPOON IMPLANT BSJO5



2087 BONE CRUSHER



2083
BONE CRUSHER
(MORSELIZER)

IMPLANT & SURGICAL | BONE INSTRUMENTS



2037T BONE GAUGE 4-6-8-10-12-14-16mm TITANIUM



2091 BONE COLLECTOR



2093 BONE COLLECTOR CHISEL STR



2093B BONE COLLECTOR STR BLADE ONLY



2093A BONE COLLECTOR CHISEL CVD



2093AB BONE COLLECTOR CVD BLADE ONLY



2092 BONE COLLECTOR WITH FILTER



2092F BONE COLLECTOR FILTER ONLY



2094
BONE MIXING BOWL



7425 BONE TRAY
SQUARE 65x58x13mm



7426 BONE TRAY BOWL
58mm DIAx28mm DEPTH

BONE SYRINGES

Bone syringes are specialized tools used in oral and maxillofacial surgery, implantology, and periodontal procedures to handle and deliver bone graft materials to specific surgical sites.

2079 BONE SYRINGE 3mm



1941 BONE SYRINGE 5mm



2080 BONE SYRINGE 3.5mm



2081 BONE SYRINGE 4mm



2082 BONE SYRINGE 6mm




MICRO SURGERY

Microsurgery instruments are precision tools used in minimally invasive procedures within dentistry, particularly in endodontics, periodontics, and implantology. These instruments enable surgeons to perform highly accurate operations on delicate oral tissues, ensuring better outcomes with minimal trauma and faster healing.



2073 MICRO SURGERY FORCEP 75mm




2074 MICRO SURGERY FORCEP 95mm




2075 MICRO SURGERY FORCEP 110mm



5546-3 MICRO SURGERY SCISSOR 18cm




5929CS SUPER CUT SURGICAL SCISSOR CASTRO CVD 14cm



5929S SUPER CUT SURGICAL SCISSOR CASTRO STR 14cm

IMPLANT & SURGICAL | MICRO SURGERY



5940 SURGICAL SCISSOR CASTRO MICRO 45° 18cm



5944 SURGICAL SCISSOR CASTRO MICRO 45° 18cm OFF ANGLE



6345 SKIN HOOK 28K-L



6346 SKIN HOOK 28K-R



6505L FORCEPS STR TEETH 1X2 15cm/6"



6506L FORCEPS CVD TEETH 1X2 15cm/6"



6511N FORCEPS MICRO SUTURE CVD 16cm



6525G FORCEPS GERALD STR MICRO VERSION 18cm



6525TC FORCEPS GERALD CVD MICRO VERSION 18cm T/C



6934 SCALPEL HANDLE BLACK EDITION FOR MICRO BLADE 15cm X 10 Dia



7730 TITANIUM, MICRO SURGICAL GERALD STR (7.1in=18cm) NON SERRATED



7731 TITANIUM, MICRO SURGICAL GERALD CVD (7.1in=18cm) NON SERRATED

ROOT FORCEPS

Root forceps are used to extract teeth and roots, and can be used on both the upper and lower gums. They are designed to prevent fractures and are effective on roots due to their serrated tips, curved jaws, and slim profile.



5105 ROOT FRAGMENTS STR



5106
ROOT
FRAGMENTS
CVD 45°



5106E
ROOT
FRAGMENTS
45°



5107
ROOT
FRAGMENTS
CVD 90°



5107E
ROOT
FRAGMENTS
90°



5108
PEET
SHARP
CVD 15°



5109
PEET
SHARP
CVD 45°



5110
PEET
SHARP
CVD 90°



5623
MICRO
ROOT
FORCEPS



4816
65
UPPER
ROOT



4817
69
UPPER
LOWER
FRAGMENT



4832
286
UPPER
ROOT



4919
30
UPPER
ROOT



4920
31
LOWER
ROOT



4922
33
LOWER
ROOT



4927
44
UPPER
ROOT



4928
45
LOWER
ROOT

IMPLANT & SURGICAL | ROOT FORCEPS



4929

51
UPPER
ROOT



4930

51A
UPPER
ROOT



4938

68
LOWER
ROOT



4964

233
LOWER
ROOT



4990

F17
LOWER
ROOT



4991

F18
UPPER
ROOT



4991T

F18
UPPER ROOT
FRAGMENT
NARROW



4992

F19
LOWER
ROOT



4992NW

F19
LOWER
ROOT
FRAGMENT
NARROW



4994

46XL
LOWER
ROOT



4994P

46XL
LOWER
ROOT



4994T

47SR
SPECIAL
ROOT
LOWER



4995

F20
UPPER
ROOT
UNIVERSAL



4996

49
UPPER
ROOT



4997

51XL
UPPER
ROOT



4997P

51XL
UPPER
ROOT



4997T

51SR
SPECIAL
ROOT
UPPER



5014

B6
LOWER
MOLAR



5015

B7
LOWER
ROOT



5016

B137
UPPER
ROOT



5018

F5
LOWER
ANTERIOR



5101

301
LOWER
ROOT



5102

300
UPPER
ROOT

TISSUE PLIERS

Manage your tissue effectively and precisely with our Tissue Pliers. Every type of tissue plier, for all types of soft tissue management.



6501 ADSON STR 12cm



6501A ADSON CVD TEETH 12cm



6501L ADSON STR TEETH 15cm



6502 ADSON STR TEETH 1x2 12cm



6502A ADSON CVD TEETH 1x2 12cm

IMPLANT & SURGICAL | TISSUE PLIERS



6502AL ADSON CVD TEETH 1x2 15cm



6502B ADSON BROWN 12cm



6502BL ADSON BROWN 15cm



6503 SEMKEN-TAYLORS STR REGULAR 12.5cm



6505 SEMKEN-TAYLORS STR TEETH 1X2 12.5cm



6505L MICROSURGICAL STR TEETH 1X2 15cm



6506 SEMKEN-TAYLORS CVD TEETH 1X2 12.5cm



6506L MICROSURGICAL CVD TEETH 1X2 15cm



6507 FORCEPS POINT 12cm

IMPLANT & SURGICAL | TISSUE PLIERS



6508 FORCEPS POINT 11.5cm



6509 FORCEPS POINT CVD 11.5cm




6510 FORCEPS POINT CVD 11.5cm



6511 SUTURE CORN TISSUE FORCEPS CVD 16cm



6511E SUTURE CORN TISSUE FORCEPS ATRAUMATIC CVD 18cm
BLACK EDITION



6511N SUTURE CORN TISSUE FORCEPS MICRO SUTURE CVD 16cm



6511S SUTURE CORN TISSUE FORCEPS STR 16cm



6512 GOLDMAN-FOX 15cm L



6513 GOLDMAN-FOX 15cm R

IMPLANT & SURGICAL | TISSUE PLIERS



6515 SELF-LOCKING 15cm



6517 LONDON-COLLEGE



6517C LONDON-COLLEGE TITANIUM



6517L LONDON-COLLEGE 15cm



6518 MERIAM



6519 PERRY 13cm




6520 MERIAM SELF-LOCKING 16cm



6523 SELF-LOCKING

IMPLANT & SURGICAL | TISSUE PLIERS



6525 GERALD STR 18cm



6525G GERALD STR MICRO VERSION 18cm



6525TC GERALD CVD 18cm MICRO



6526 TOMOR 16cm



6526L TOMOR 20cm



6527 TP16



6528 TP50L



6531 DEBAKEY ATRAUMATIC STR 18cm



6531B DEBAKEY ATRAUMATIC STR 18cm BLACK TITANIUM



6531C DEBAKEY ATRAUMATIC CVD 18cm

IMPLANT & SURGICAL | TISSUE PLIERS



6531CB DEBAKEY ATRAUMATIC CVD 18cm BLACK TITANIUM



6532 GERALD STR REGULAR 18cm



6532C GERALD CVD REGULAR 17.5cm



6533 SPECIAL MEMBRANE GRAFTING



6534 ABUTMENT HOLDER
HOLE SIZE 1-2.5mm/1.6-3.2mm



7730 TITANIUM, MICRO SURGICAL GERALD STR (7.1in=18cm) NON SERRATED



7731 TITANIUM, MICRO SURGICAL GERALD CVD (7.1in=18cm) NON SERRATED



6442 CHALAZION 9mm



6443S CHALAZION 16x10.5mm



6443 CHALAZION 24x14mm



6444 CHALAZION 24x14mm

CALIPERS & MEASURING

Calipers are used to take accurate measurements in specific places around the bone. They are also used for measuring the size and shape of teeth to ensure accurate restorations, such as crowns, bridges, or dentures.



4604
IVANSON
CALIPER



4606
BOLEY
CALIPER 14cm



4613
CASTRO
CALIPER
8.5cm




4613C
CASTRO
CALIPER
CVD 8.5cm



4613L
CASTRO
CALIPER CLC
40L STR 17cm



4629
FINGER RULER 38mm



4618 RULER 15cm 6in

IMPLANT & SURGICAL | CALIPERS & MEASURING



4613LA
CASTRO
CALIPER CLC
40C CVD 17cm



4613RA
CASTRO
CALIPER WEISS
MODIFIED CLC
CVD 40 CWM 17cm




4614
BONE RIDGE
MAPPING
CALIPER
0-35mm



4625
BONE RIDGE
MAPPING
IMPLANT
CALIPER




4628 MICRO BOLEY CALIPER



4650 IMPLANT CALIPER



4654 VISAGE MEASUREMENT CALIPER



4655 ZYGOMA CALIPER

TISSUE PUNCHES

Remove soft tissue with precision and improve accuracy and access in implant placement with our Dental Tissue Punches.



2076 STR 4mm



2077 STR 5mm



2078 STR 6mm



1955A CVD 2mm



1956A CVD 3mm



2076A CVD 4mm



2077A CVD 5mm



2078A CVD 6mm

ACCU BERNARD SPADE SET

Our Accu Bernard Spade Set is manufactured with surgical grade stainless steel coated with blue titanium. The fine micro serrations easily loosen the tooth from the socket and periodontal ligaments. The syndesmotomes can be used to elevate small roots.




3771 ACCU LUX-TOME SPADE



3790 ACCU BERNARD SPADE OFFSET IN



3791 ACCU BERNARD SPADE OFFSET OUT



3792 ACCU BERNARD SPADE LEFT



3793 ACCU BERNARD SPADE RIGHT

PERIOSTEALS

Periosteals are used to separate and elevate the periosteum after removing gingival tissue. They can also be used for tooth extraction, bone grafting and sinus lifts.



3501 BENNETT



3502 P149



3506 PERIOSTEALS MOLT9



3503 HOPKINS1



3505M MEAD3 SMALL



3504 HOPKINS2



3506E PERIOSTEALS MOLT9 BLACK EDITION



3505 MEAD3



3506T PERIOSTEALS MOLT9 TITANIUM

IMPLANT & SURGICAL | PERIOSTEALS



3506M PERIOSTEALS MOLT9 MINI




3507 GF14



3516 BUSER



3508 3/9A



3516E BUSER BLACK EDITION



3509 PR3



3516SP BUSER SP



3509T PR3 TITANIUM



3516T BUSER TITANIUM



3510 P16



3517 KN1



3511 P15



3518 KN152



3512 HOWARD 9



3521 HF20



3512S HOWARD 9 MINI



3521T HF20 TITANIUM




3513 ALLEN9



3524 FREER



3513B ALLEN9 BLACK EDITION



3526 WOODSON1 MODIFY



3513T ALLEN9 TITANIUM



3526T WOODSON1 MODIFY TITANIUM



3514 24G



3534T PERIO SURGICAL TITANIUM



3514T 24G TITANIUM



3536 P1



3515 BENQUE



3537 P2

IMPLANT & SURGICAL | PERIOSTEALS



3520 HOURIGAN2



3522 MOLT5



3523 SELDIN23



3523N SELDIN23 BLUNT



3525 OHL10

IMPLANT & SURGICAL | PERIOSTEALS



3531 P25



3531N P25 NEW HANDLE



3532 P25L



3533 P25XL



3533E P25XL BLACK EDITION



3533N P25XL 6EZ



3540 PRICHARD DOUBLE END TITANIUM

SURGICAL CURETTES (SPOON)

Surgical curettes are used for curettage and cyst removal as well as craping or debriding tissue. The spoon tip allows for precision in removing unwanted or unhealthy tissue without damaging surrounding areas.



4101 LUCAS84



4101S LUCAS84S SERRATED



4103AS LUCAS86AS SERRATED
OFF ANGLE



4101SE LUCAS84S SERRATED
BLACK EDITION



4103ASE LUCAS86AS
SERRATED OFF ANGLE
BLACK EDITION



4102 LUCAS85



4103S LUCAS86S SERRATED



4102S LUCAS85 SERRATED



4103SE LUCAS86S SERRATED
BLACK EDITION



4102SE LUCAS85S SERRATED
BLACK EDITION



4107 LUCAS87




4103 LUCAS86



4108 LUCAS88

IMPLANT & SURGICAL | SURGICAL CURETTES (SPOON)



4104 MILLER8



4105 MILLER9



4106 MILLER10



4109 MILLER11



4110 MILLER12



4111 MOLT2/4



4126

POWER SERRATED CURETTE KIT



4101S LUCAS84S SERRATED



4102S LUCAS85 SERRATED



4103AS LUCAS86AS SERRATED OFF ANGLE



4103S LUCAS86S SERRATED

4126E

POWER SERRATED CURETTE KIT BLACK EDITION



4101SE LUCAS84S SERRATED BLACK EDITION



4102SE LUCAS85S SERRATED BLACK EDITION



4103SE LUCAS86S SERRATED BLACK EDITION



4103ASE LUCAS86AS SERRATED OFF ANGLE BLACK EDITION

BONE CHISELS

Bone chisels are used to shape and smooth bone tissue around teeth. They are also used to remove excess bone, trim alveolar bone, split bone, and prepare the bone for dental implants.



4201 GAR1



4202 GAR2




4209 C4




4210 W1




4211 W2



4212 M18



4213 C7



4214 C8



4215 C9

BONE FILES

Bone files are used to remove small bone segments, tooth decay, and shape dental fillings. They can also be used to smooth bone and narrow the canals of the interior of a tooth. Their serrated edges make it easy to maneuver between soft tissue and bone.



4216 45



4217 1W



4218 2W



4219 1X



4220 2X



4221 3X



4222 21

MALLETS

Mallets are used to impact dental materials and instruments during dental procedures. They are also used to split teeth, remove bone, and break up bone tissue. They are often used with chisels or wedges.



4223N MEAD1



4223M MEAD MINI


4223MT MEAD MINI REPLACEMENT TIPS (2)
4223T MEAD REPLACEMENT TIPS (2)




4226 MEAD MINI
FOR IMPLANT

CHISELS & GOUGES

Chisels are used to remove, smooth, and contour bone so that the implant can be placed. Gouges are used to loosen teeth from the gums to make extraction easier and to prevent damage to surrounding tissue and teeth.



4503 FREER 16cm 4mm CVD



4504 LUCAS 16cm 4mm



4507 14cm 6mm



4508 PARTCSH 14cm 5mm



4510 PARTCSH 14cm 4mm

BONE RONGEURS

Bone rongeurs are used to cut or remove small pieces of bone or tissue. In addition to being used to trim and contour bone they can also be used to open a window in the bone to access tissue underneath.



4401
BLUMENTHAL
15cm



4402
BLUMENTHAL
45° 14.5cm



4403
BLUMENTHAL
90° 14.5cm



4405
BEYER
DOUBLE
ACTION
18cm



4406
BEYER
DOUBLE
ACTION
RUSKIN-MINI
15.5cm



4408
FRIEDMAN
14cm



4409
MINI
FRIEDMAN 15
12cm



4411
MICRO
FRIEDMAN 15
14.5cm



4411L
MICRO
FRIEDMAN 15
15.5cm



4412
MICRO
FRIEDMAN 90
14.5cm



4413
MICRO
FRIEDMAN 45
14.5cm D-42

HEMOSTATIC FORCEPS

Hemostatic forceps are used to clamp blood vessels, control bleeding, and move gauze and other materials around the mouth. They can also be used to hold and move small objects like bone and tooth fragments.



5301C
BABY
MOSQUITO
CVD 10cm



5302C
MICRO-
HALSTED
CVD 12.5cm



5303C
HALSTED-
MOSQUITO
CVD 12.5cm

5301
BABY MOSQUITO STR 10cm

5302
MICRO-HALSTED STR 12.5cm

5303
HALSTED-MOSQUITO STR 12.5cm

IMPLANT & SURGICAL | HEMOSTATIC FORCEPS



5305C
MICRO-
HALSTED
CVD 12.5cm

5305
HALSTED-MOSQUITO STR 12.5cm



5307
PEAN 12.5cm



5309C
PEAN
BABY-CRILE
CVD 14cm H-2

5309
PEAN BABY-CRILE STR 14cm H-1



5310C
CRILE/
KELLY
CVD 14cm

5310
CRILE/KELLY STR 14cm



5311C
BABY CRILE/
KELLY KOCHER
CVD 14cm

5311
BABY CRILE/KELLY KOCHER
STR 14cm



5312C
CHILE
KOCHER
CVD 16cm

5312
CHILE KOCHER STR 16cm



5313C
HALSTED
CVD 18.5cm

5313
HALSTED STR 18.5cm



5315C
HALSTED
CVD 18.5cm

5315
HALSTED STR 18.5cm



5316
BABY-MIXTER CVD 18cm



5318
IMPLANT ABUTMENT
HOLDER CVD 14.5cm

MATHIEU PLIERS

Mathieu pliers are used to hold, place, and tie elastics and ligatures, and to guide needles during suturing, are lightweight and flexible, with a spring function that allows the user to lock and unlock the jaws.



5407 MATHIEU
14cm



5407F MATHIEU
FINE 14cm



5407U MATHIEU
US ORTHO 14cm



5414 MATHIEU
DELICATE STR 15cm



5415 MATHIEU
DELICATE CVD 15cm



5418 MATHIEU
14cm



5419 MATHIEU
THIN 14cm



5501N MATHIEU T/C THIN
BLACK TITANIUM 16cm



5502 MATHIEU
T/C WIDE 14cm



5502L MATHIEU
T/C WIDE 17cm

NEEDLE HOLDERS

Needle holders are used to grasp and manipulate needles during suturing. Needle holders can also be used in other dental procedures where precise needle control is required.



5402 CRILE-MURRAY 15cm



5403
MAYO-
HEGAR
15cm



5403L
MAYO-
HEGAR
16cm



5403W
MAYO-
HEGAR
WIDE
15cm



5408 BOYNTON 12cm



5409B CASTROVIEJO STR
BLACK TITANIUM 16cm



5409C CASTROVIEJO CVD T/C 16cm

IMPLANT & SURGICAL | NEEDLE HOLDERS



5505L

CRILE-WOOD T/C 18cm



5505TN

CRILE-WOOD STR T/C 15cm



5506

CRILE-WOOD CVD T/C 15cm



5508

RYDER STR T/C 18cm



5509

DERF STR T/C 11.5cm



5512

HALSEY STR T/C 13cm



5521

D-47 MAYO-HEGAR STR T/C 18cm



5521L

MAYO-HEGAR STR T/C 20cm



5545

CRILE-WOOD "S" CVD HANDLES T/C 15cm

IMPLANT & SURGICAL | NEEDLE HOLDERS



5504B OLSEN-HEGAR
BLACK TITANIUM STR 15cm



5505B CRILE-WOOD
BLACK EDITION STR 15cm



5505BN CRILE-WOOD
BLACK TITANIUM STR 15cm



5505LB CRILE-WOOD
BLACK EDITION STR T/C 18cm



5505TB CRILE-WOOD
BLACK TITANIUM
STR T/C 15cm



5505TNB CRILE-WOOD
BLACK EDITION
STR T/C 15cm



5506B CRILE-WOOD
BLACK EDITION
CVD T/C 15cm



5509B DERF
BLACK EDITION
STR T/C 11.5cm



5512B HALSEY
BLACK EDITION
STR T/C 13cm



5520B MAYO-HEGAR
BLACK EDITION
STR T/C 16cm



5521B MAYO-HEGAR
BLACK TITANIUM
STR T/C 18cm

IMPLANT & SURGICAL | NEEDLE HOLDERS



5513B

CASTROVIEJO
BLACK EDITION
STR 16cm



5513BN

CASTROVIEJO
BLACK TITANIUM
STR 16cm



5514B

CASTROVIEJO
BLACK EDITION
STR 18cm



5514BN

CASTROVIEJO
BLACK TITANIUM
STR 18cm



5515B

CASTROVIEJO
BLACK EDITION
CVD 16cm



5515BN

CASTROVIEJO
BLACK TITANIUM
CVD 16cm



5516B

CASTROVIEJO
BLACK EDITION
CVD 18cm



5516BN

CASTROVIEJO
BLACK TITANIUM
CVD 18cm



5527

DOUBLE
ACTION JAWS
14cm



5531

CASTROVIEJO
SPECIAL
CVD 45°
16.5cm

SCISSORS

Scissors are used to cut tissue, sutures and soft material like gauze and membranes.



5905
KELLY
CVD
16cm

5904 KELLY STR 16cm



5906C
QUINBY
CVD
12.5cm

5906 QUINBY STR 12.5cm



5907 DEAN 17.5cm



5907N DEAN
GOLD TITANIUM
SCREW 17.5cm

IMPLANT & SURGICAL | SCISSORS



5908C
METZENBAUM
CVD
14.5cm

5908 METZENBAUM
STR 14.5cm



5909 LOCKLIN 16cm



5911C
GOLDMAN-FOX
CVD
12.5cm

5911 GOLDMAN-FOX
STR 12.5cm



5916 GOLDMAN-FOX
SPECIAL CVD 12.5cm



5917 LA GRANGE
CVD 11.5cm



5919 IRIS STR 11.5cm



5920
IRIS
CVD
11.5cm




5920F IRIS MICRO
CVD 11.5cm



5922 SUTURE 15cm

IMPLANT & SURGICAL | SCISSORS



5905B
KELLY
SUPER CUT
BLACK EDITION
CVD 16cm



5907B DEAN SUPER CUT
BLACK EDITION 17.5cm



5907BN DEAN
POWER BOW 17.5cm



5908CB
METZENBAUM
BLACK EDITION
CVD 14cm

5908B METZENBAUM
BLACK EDITION STR 14.5cm



5911CB
GOLDMAN-FOX
BLACK EDITION
CVD 12.5cm

5911B GOLDMAN-FOX
BLACK EDITION STR 12.5cm



5911CBN
GOLDMAN-FOX
POWER BOW
CVD 12.5cm

5911BN GOLDMAN-FOX
POWER BOW STR 12.5cm



5916B GOLDMAN-FOX
SPECIAL CVD
BLACK EDITION 12.5cm



5917B LA GRANGE
SUPER CUT
BLACK EDITION CVD 11.5cm



5917BN LA GRANGE
POWER BOW CVD 11.5cm




5920BN
IRIS
POWER BOW
CVD 11.5cm

5919BN IRIS
POWER BOW STR 11.5cm



5919B IRIS SUPER CUT
BLACK EDITION STR 11.5cm



5920B
IRIS
SUPER CUT
BLACK EDITION
CVD 11.5cm



5922B SUTURE SUPER CUT
BLACK EDITION 15cm D-52

IMPLANT & SURGICAL | SCISSORS



5905S
KELLY
SUPER CUT
CVD 16cm

5904S KELLY SUPER CUT
STR 16cm



5907S DEAN
SUPER CUT 17.5cm



5908SC
METZENBAUM
SUPER CUT
CVD 14.5cm

5908S METZENBAUM
SUPER CUT STR 14.5cm



5908SCP
METZENBAUM
POINT FRONT
CVD 14.5cm

5908SP METZENBAUM
POINT FRONT STR 14.5cm



5911S
GOLDMAN-FOX
SUPER CUT
CVD 12.5cm

5911S GOLDMAN-FOX
SUPER CUT STR 12.5cm



5916S GOLDMAN-FOX
SUPER CUT TUNGSTEN
BLACK EDITION 12.5cm



5920S
IRIS
SUPER CUT
CVD 11.5cm

5919S IRIS
SUPER CUT STR 11.5cm

IMPLANT & SURGICAL | SCISSORS



5923
SUTURE SLIM
12cm



5924
CROWN
STR 12cm



5926
CASTRO
STR 10cm



5926B
CASTRO
BLACK
TITANIUM
STR 10cm



5926C
CASTRO
CVD 10cm



5929
CASTRO
STR 14cm



5929C
CASTRO
CVD 14cm



5929CS
CASTRO
MICRO
CVD 14cm



5929S
CASTRO
MICRO
STR 14cm



5940
CASTRO
MICRO
LOCKLIN
45° 18cm



5944
CASTRO
MICRO
OFF ANGLE
45° 18cm

SUTURE SCISSORS

Suture scissors feature a hook design specifically utilized to lift the suture and cut with a single sweep.



5904B KELLY SUPER CUT
BLACK EDITION STR 16cm



5905S KELLY SUPER CUT
CVD 16cm



5907B DEAN SUPER CUT
BLACK EDITION 17.5cm



5907S DEAN
SUPER CUT 17.5cm

IMPLANT & SURGICAL | SUTURE SCISSORS



5908SC METZENBAUM
SUPER CUT CVD 14.5cm



5911B GOLDMAN-FOX
BLACK EDITION STR 12.5cm



5911CB GOLDMAN-FOX
BLACK EDITION CVD 12.5cm



5911S GOLDMAN-FOX
SUPER CUT STR 12.5cm



5911SC GOLDMAN-FOX
SUPER CUT CVD 12.5cm



5917B LA GRANGE
SUPER CUT
BLACK EDITION CVD 11.5cm



5919S IRIS
SUPER CUT STR 11.5cm



5922 SUTURE 15cm



5922B SUTURE SUPER CUT
BLACK EDITION 15cm



5923
SUTURE SLIM
12cm



5923A
SUTURE SLIM
CVD 12cm



5944
CASTRO
MICRO
OFF ANGLE
45° 18cm

TUNGSTEN CARBIDE INSERT SCISSORS

Harder than traditional stainless steel, cut through the toughest of materials with our Tungsten Carbide Scissors.



6001 METZENBAUM
ROUND T/C STR 14.5cm



6002
METZENBAUM
ROUND
T/C CVD 14.5cm



6003 METZENBAUM
POINT T/C STR 14.5cm



6004
METZENBAUM
POINT
T/C CVD 14.5cm



6003BN METZENBAUM
POWER BOW
T/C STR 15cm



6004BN
METZENBAUM
POWER BOW
T/C CVD 14.5cm



6008 WIRE CUTTING
T/C 12cm



6009 IRIS STR
T/C 11.5cm



6010
IRIS CVD
T/C 11.5cm



6011 DEAN
T/C 17.5cm



6012 LOCKLIN
T/C 16.5cm



6014 KELLEY STR
T/C 16cm



6015
KELLEY CVD
T/C 16cm

SKIN HOOKS


Skin hooks are used hold back the edges of skin during dental procedures. They can also be used to elevate tissue and mobilize it as well as provide gentle horizontal traction.



6312 GILLIES 18cm



6345 MICRO SURGERY 28K-L



6346 MICRO SURGERY 28K-R

RETRACTORS

Retractors are used to hold back soft tissues, such as gums and cheeks, to provide better visibility and access to the treatment area. Mouth gags are used to prop the mouth open during procedures.

6342 FARABEUF 12.5cm



6343 FARABEUF 15cm



6344 ROUX 16.5cm



6420A UNIVERSITY OF MINNESOTA
FLAT WITH SERRATION



6420 UNIVERSITY OF MINNESOTA



6420B UNIVERSITY OF MINNESOTA
CONCAVE FRONT



6420N UNIVERSITY OF MINNESOTA
SPECIAL HANDLE HEAVY HOLD



6420R UNIVERSITY OF MINNESOTA
ROUND FRONT



6421 ALM 7cm



IMPLANT & SURGICAL | RETRACTORS



6423 CAWOOD MINNESOTA



6428 BISHOP



6426 MOUTH GAG ADULT

6426T BLUE REPLACEMENT TIPS (ADULT SIZE)


6427 MOUTH GAG CHILD

6427S MOUTH GAG PEDO

6427T BLUE REPLACEMENT TIPS (CHILD SIZE)



6424 COLUMBIA 14cm



6425 ORRINGER HANDS FREE




6429 HUMAN ADULT



6430 HUMAN CHILD



6433 VESTIBULUM LIP LARGE



6434 VESTIBULUM LIP SMALL



6434N UNIVERSAL CHEEK
AND LIP SMALL



6435 BAYONET 15cm

IMPLANT & SURGICAL | RETRACTORS



6438 SPANDEX
SMALL 10cm



6440 S/E 9cm



6446
MOUTH PROPS
ADULT



6447
MOUTH PROPS
CHILD



6445 LIP RETRACTOR



6436
WEIDER TONGUE
RETRACTOR 30mm



6436L
WEIDER TONGUE
RETRACTOR 35mm



6449 STERNBERG RETRACTOR



6448 AUSTIN



6448A AUSTIN DOUBLE



6450 CHEEK RETRACTOR



6451 CHEEK DOUBLE SIDE



6452 SINUS SOFT TISSUE HOLD AND EXPANDER 10mm



6453 SINUS SOFT TISSUE HOLD AND EXPANDER 15mm



6454 CHEEK AND TONGUE SPOON TYPE

SCALPEL HANDLES

Scalpel handles are used in conjunction with mess blades for making precise incisions in soft tissue. Our handles are designed for comfort and control during procedures.



6901 4 FLAT




6902 3G FLAT



6904 5A ROUND CVD



6905 5A ROUND STR



6905D 5A ROUND DOUBLE HEAD



6906 7 THIN HANDLE

IMPLANT & SURGICAL | SCALPEL HANDLES



6907 3L FLAT LONG




6920 6EZ



6920A PALATAL HARVESTING 1.5mm GAP/1mm



6920B PALATAL HARVESTING 1.5mm GAP/2mm




6920BE IMPLANT SURGERY BLACK EDITION



6920G 6EZ DR. GEORGE TIP SHAPE 45° CUT



6924 S/E 360° BLADE ANGLE TURN AROUND



6934 BLACK EDITION FOR MICRO BLADE



6935 D/E 45° ORBAN TYPE OFF ANGLE
S/E 360° BLADE ANGLE TURN AROUND

IMPLANT & SURGICAL | SCALPEL HANDLES


6920S PALATAL HARVESTING KIT




6920A SCALPEL HANDLE



6920B SCALPEL HANDLE



6966 EASY BLADE REMOVAL FORCEP




6 way adjustable angled
Scalpel Handle



6923
SCALPEL HANDLE
180° ADJUSTABLE
16.2cm

SYRINGES

Syringes are used for injecting anesthetics from a cartridge. Three-way syringes supply air, water or a mist through separate interval channels created by combining pressurized air with water flow.



7001 INTRALIGMENTARY 1.8ml



7003 1.8ml WITH HARPOON



7017 PAROJECT

IMPLANT & SURGICAL | SYRINGES



7018L ASPIRATED 2.2ml WING



7018LN ASPIRATED 2.2ml WING



7018N ASPIRATED 1.8ml WING




7019 3-WAY



7019T 3-WAY TIPS (2)



ELEVATORS & LUXATORS




E

ELEVATORS

Elevators are used in a variety of dental procedures, including routine extractions, removing impacted teeth, and loosening teeth that are malpositioned or at risk of fracture. They work on the principle of leverage to dislodge a tooth from its socket.

(SHOWN ACTUAL SIZE)



3701 HEIDBRINK 1



3702

HEIDBRINK 2



3704

SELDIN 12M



3705

SELDIN 34



3706

SELDIN 34S



3706S

SELDIN 34S
SERRATED



3707

HEIDBRINK 40



3708

HEIDBRINK 41



3709

APEXO 80

ELEVATORS & LUXATORS | ELEVATORS



3710
APEXO 81



3711
APEXO 81E



3712
APEXO 301



3712S
APEXO 301
SERRATED



3712W
APEXO 301
WIDE



3712WS
APEXO 301
WIDE
SERRATED



3713
APEXO 302
CVD INNER 2.8mm
±5%



3714
SELDIN 2S



3715
APEXO 303
CVD INNER 2.8mm
±5%



3716
CRANE PICK
8



3717
SELDIN 1L



3718
SELDIN 1R



3719
SELDIN 4L



3720
SELDIN 4R



3721
CRYER M25



3722
CRYER M26



3725
CRYER S30A



3726
CRYER S31A



3729
CRYER L46L



3730
CRYER L46R



3731
BARRY 50



3732
BARRY 51



3734
MILLER
APEXO 71



3735
MILLER
APEXO 72

ELEVATORS & LUXATORS | ELEVATORS

SHOWN ACTUAL SIZE



3745 FLOHR R-SHARP

SHOWN SMALLER THAN ACTUAL SIZE



3747S BERNARD 3 LEFT SMALL HANDLE



3750S BERNARD 4 RIGHT SMALL HANDLE



3736
MILLER
APEXO 73



3737
MILLER
APEXO 70



3738
MILLER
APEXO 74



3739
APICAL 2S/301
EXTRA THIN



3740
APICAL 2S/302
EXTRA THIN



3741
APICAL 4R/303
EXTRA THIN



3742
FLOHR
L-SHARP



3743
FLOHR
L-SHARP



3747
BERNARD 3
LEFT



3748
BERNARD 2
STR



3749
BERNARD 2
CVD



3749IN
BERNARD
IN



3749OUT
BERNARD
OUT



3750
BERNARD 4
RIGHT



3750D
FLAP
OPENER



3750DT
FLAP
OPENER
TITANIUM GOLD

ELEVATORS & LUXATORS | ELEVATORS



3759
COPELAND 3
(4.5mm)



3760
COPELAND 2
(3.5mm)



3765
LINDO
STR 3mm



3766
LINDO
STR 4mm



3767S
SELDIN
46S IN
SERRATED



3768S
SELDIN
46S OUT
SERRATED



3787
FRIEDMAN
GOUGE
F31 OUT
SERRATED



3788
FRIEDMAN
GOUGE
F31 IN
SERRATED



3801
FLOHR
3R



3802
FLOHR
1S



3803
FLOHR
2L



3805
KOPP
1S



3806
KOPP
3R



3807
WWJL
SMALL HANDLE



3808
WWJS
SMALL HANDLE



3809
WWJR
SMALL HANDLE



3814
CRYER
39



3814S
CRYER 39
SMALL HANDLE



3815
CRYER
40R



3815S
CRYER 40R
SMALL HANDLE



3818S
E77R OUT
4mm
SERRATION



3820S
E77 INNER
4mm
SERRATION



3822
HYLIN
2S



3823
HYLIN
1S



3824
HYLIN
3S



3827 COLEMAN 2LS



3828
COLEMAN
2RS



3830
LINDO LEVIN
INNER
SERRATED



3831
LINDO LEVIN
OUTER
SERRATED



3841
COGSWELL
B



3842
WING 1
EXTRA THIN



3843
WING 2
EXTRA THIN



3844
WING
3



3845
WING
4



3852
POTTS
6



3853
POTTS
7

ELEVATORS & LUXATORS | ELEVATORS



3837 BARRY 17L T-BAR



3838 BARRY 17R T-BAR



3839 WINTER 14L T-BAR



3840 WINTER 14R T-BAR



3850 HOURIGAN 4 T-BAR



3851 HOURIGAN 3 T-BAR



3852T POTTS 6 T-BAR



3853T POTTS 7 T-BAR

ELEVATORS & LUXATORS | ELEVATORS



3856
ASL
2.5mm L



3857
ASL
2.5mm R



3878
HYLIN
12cm
EXTRA THIN



3878L
HYLIN
16cm
EXTRA THIN



3882 LIZA L



3883 LIZA R

ELEVATORS & LUXATORS | ELEVATORS

4034

POWER FLEX ELEVATOR SET



4034-1 POWER FLEX APEXO 301 STR 2.8mm ±5% TITANIUM



4034-2 POWER FLEX APEXO 302 CVD IN 2.8mm ±5% TITANIUM



4034-3 POWER FLEX APEXO 303 CVD OUT 2.8mm ±5% TITANIUM



4034-4 POWER FLEX SELDIN 34S 4.0mm ±5% TITANIUM



4034-5 POWER FLEX E77 UP 4.0mm ±5% TITANIUM



4034-6 POWER FLEX E77 DOWN 4.0mm ±5% TITANIUM

4031

POWER TWIST ELEVATOR SET



4023 POWER TWIST ELEVATOR STR-L



4024 POWER TWIST ELEVATOR STR-R



4026 POWER TWIST ELEVATOR CVD-L



4027 POWER TWIST ELEVATOR CVD-R

LUXATORS

Luxators are used to loosen and remove teeth from the mouth. They are sharp with a thin blade that's inserted into the periodontal ligament to cut the fibers and separate the tooth from the socket.



3751 2mm STR



3752
3mm STR



3753
5mm STR



3754
2mm CVD



3755
3mm CVD



3756
5mm CVD



3776
LUX-TOME
3mm



3777
LUX-TOME
4mm



3778
LUX-TOME
5mm



3779
3mm STR
SMALL HANDLE



3780
4mm STR
SMALL HANDLE



3781
5mm STR
SMALL HANDLE



3782
3mm CVD
SMALL HANDLE



3783
4mm CVD
SMALL HANDLE



3784
5mm CVD
SMALL HANDLE



3785
4mm STR
SMALL HANDLE



3786
4mm CVD
SMALL HANDLE

ELEVATORS & LUXATORS | LUXATORS

3799

POWER FLEX PDL TITANIUM SET OF 7



3799-1 POWER FLEX PDL TITANIUM ELEVATOR 3mm STR



3799-2 POWER FLEX PDL TITANIUM ELEVATOR 3mm CVD



3799-3 POWER FLEX PDL TITANIUM ELEVATOR SPADE



3799-4 POWER FLEX PDL TITANIUM SMALL DISTAL



3799-5 POWER FLEX PDL TITANIUM SMALL MESIAL



3799-6 POWER FLEX PDL TITANIUM LARGE DISTAL



3799-7 POWER FLEX PDL TITANIUM LARGE MESIAL

ELEVATORS & LUXATORS | LUXATORS

4035


POWER FLEX LUXATOR SET



4035-1 POWER FLEX LUXATOR STR 2.5mm ±5% TITANIUM



4035-2 POWER FLEX LUXATOR CVD 2.5mm ±5% TITANIUM



4035-3 POWER FLEX LUXATOR STR 3.5mm ±5% TITANIUM



4035-4 POWER FLEX LUXATOR CVD 3.5mm ±5% TITANIUM



4035-5 POWER FLEX LUXATOR STR 4.5mm ±5% TITANIUM



4035-6 POWER FLEX LUXATOR CVD 4.5mm ±5% TITANIUM

LUX-TOMES

When needing the ability to create space and luxate, but also have enough thinness to insert deep into the periodontal pockets, our Power Lux-Tomes offer a hybrid approach to performing extractions.



3769 ACCU LUX-TOME STR SERRATED TITANIUM



3770 ACCU LUX-TOME CVD SERRATED TITANIUM



3771 ACCU LUX-TOME SPADE SERRATED TITANIUM



3772 ACCU LUX-TOME DISTAL SERRATED TITANIUM



3773 ACCU LUX-TOME MESIAL SERRATED TITANIUM



3774 ACCU LUX-TOME LARGE DISTAL SERRATED TITANIUM



3775 ACCU LUX-TOME LARGE MESIAL SERRATED TITANIUM


ELEVATORS & LUXATORS | LUX-TOMES




3790 ACCU BERNARD SPADE OFFSET IN



3791 ACCU BERNARD SPADE OFFSET OUT



3792 ACCU BERNARD SPADE LEFT



3793 ACCU BERNARD SPADE RIGHT



3794 ACCU LUX-TOME T-BAR

POWER FLEX LUX-TOMES

The Power Lux-Tomes that you know and love, except enhanced with NiTi tips for flexibility, durability, and a lightweight handle design that allows for maximum control and usage inside the mouth.

3796 POWER FLEX SPADE KIT



3796-1 POWER FLEX SPADE STR



3796-2 POWER FLEX SPADE CVD



3796-3 POWER FLEX SPADE LEFT



3796-4 POWER FLEX SPADE RIGHT

ELEVATORS & LUXATORS | POWER FLEX LUX-TOMES



3797-1 POWER FLEX ACCU LUX-TOME 3mm STR SERRATED



3797-2 POWER FLEX ACCU LUX-TOME 3mm CVD SERRATED



3797-4 POWER FLEX ACCU LUX-TOME SMALL DISTAL SERRATED



3797-5 POWER FLEX ACCU LUX-TOME SMALL MESIAL SERRATED



3797-6 POWER FLEX ACCU LUX-TOME MEDIUM DISTAL SERRATED



3797-7 POWER FLEX ACCU LUX-TOME MEDIUM MESIAL SERRATED

3798

POWER FLEX PDL ACCU LUXATING FULL SET OF 10

3796-1 POWER FLEX SPADE STR

3796-2 POWER FLEX SPADE CVD

3796-3 POWER FLEX SPADE LEFT

3796-4 POWER FLEX SPADE RIGHT

3797-1 POWER FLEX ACCU LUX-TOME 3mm STR SERRATED

3797-2 POWER FLEX ACCU LUX-TOME 3mm CVD SERRATED

3797-4 POWER FLEX ACCU LUX-TOME SMALL DISTAL SERRATED

3797-5 POWER FLEX ACCU LUX-TOME SMALL MESIAL SERRATED

3797-6 POWER FLEX ACCU LUX-TOME MEDIUM DISTAL SERRATED

3797-7 POWER FLEX ACCU LUX-TOME MEDIUM MESIAL SERRATED

ROOT TIP PICKS

Root tip picks are used to retrieve root tips that are in difficult to reach areas or that have broken away from the tooth during extraction. The sharp tip is placed into the periodontal ligament space and the root is teased from the socket.



3901 HEIDBRINK 1



3902

HEIDBRINK 2



3903

HEIDBRINK 3



3904

HEIDBRINK 4



3905

HEIDBRINK 5



3907

APICAL 9



3908

APICAL 9R



3916

HOWARD 1



3917

E-17 HOWARD 2



3918

HOWARD 3



3919

HOWARD 4



3920

HOWARD 5



4020

78 STR



4021

79L



4022

80R

ELEVATORS & LUXATORS | ROOT TIP PICKS



3908N APICAL 9R 6EZ



3909E APICAL 9L 6EZ



3906 HEIDBRINK 13/14



3906S HEIDBRINK 13/14 SERRATED



3921 HOWARD 2/3

4001S POWER ROOT TIP PICKS SET OF 2



4001 POWER ROOT TIP PICK LEFT DAVIS TEASER 2



4002 POWER ROOT TIP PICK RIGHT DAVIS TEASER 3

4033M POWER FLEX PREMIUM TWIST ROOT PICKS SET OF 3



4033M-1 POWER FLEX TWIST PERIOTOME/ROOT PICK ANTERIOR LEFT S/E



4033M-2 POWER FLEX TWIST PERIOTOME/ROOT PICK ANTERIOR RIGHT S/E



4033M-3 POWER FLEX TWIST PERIOTOME/ROOT PICK POSTERIOR D/E



FORCEPS



F

FORCEPS | ATRAUMATIC EXTRACTION FORCEPS

NO FRACTURE



4974
F1
UPPER
PREMOLARS

NO TISSUE DAMAGE



4975
F2
LOWER
PREMOLARS

EXTRACTION WITH EASE



4976
F2
LOWER
PREMOLARS



4977N
F4N
LOWER
ANTERIOR



4979
F6
UPPER
MOLARS



4980C
F7C
LOWER
MOLARS

AMERICAN PATTERN

American pattern dental extraction forceps are typically identified by the shape of connecting point where two forceps beaks connect. They feature hinges that are completely polished out and smooth so that the screw isn't visible.



FORCEPS | EXTRACTION FORCEPS | AMERICAN PATTERN



FORCEPS | EXTRACTION FORCEPS | AMERICAN PATTERN



ENGLISH PATTERN

English pattern forceps have visible screws in the hinge area. The hinge is vertical, the beaks are at a right angle to the handle and the fulcrum is closer to the working end.



4901 1 UPPER CENTRAL



4902
2
UPPER
BICUSPID



4905
7
UPPER
PREMOLAR



4906
8
LOWER
PREMOLAR



4908
17
UPPER RIGHT
MOLAR



4909
18
UPPER LEFT
MOLAR

FORCEPS | EXTRACTION FORCEPS | ENGLISH PATTERN




FORCEPS | EXTRACTION FORCEPS | ENGLISH PATTERN



EXTRACTION FORCEPS


Extraction forceps are primarily used for the removal of teeth that are decayed, fractured, or impacted. They are also used to remove remaining roots after a tooth has been fractured or previously extracted.



4972 22 L-L/U-R SPECIAL MOLAR



4973
22
L-R/U-L
SPECIAL
MOLAR



4974
F1
UPPER
UNIVERSAL



4974B
PEDO F1
UPPER
UNIVERSAL



4974P
F1P
UPPER
UNIVERSAL




4975
F2
LOWER
UNIVERSAL

FORCEPS | EXTRACTION FORCEPS



FORCEPS | EXTRACTION FORCEPS



FORCEPS | EXTRACTION FORCEPS



5013
B5
LOWER
INCISOR



5014
B6
LOWER
MOLAR



5015
B7
LOWER
ROOT



5016
B137
UPPER
ROOT



5017
B139
UPPER
PREMOLAR



5018
F5
LOWER
ANTERIOR



5019
PEDO F5
LOWER
ANTERIOR



5019SC
F6
UPPER
MOLARS



5020
PEDO F6
UPPER
MOLAR



5020SC
PEDO F1
UPPER
UNIVERSAL



5022
B23
LOWER
MOLAR



5101
301
LOWER
ROOT



5102
300
UPPER
ROOT



5103 EXTRACTION TWEEZER ALBRECHT 16cm ROUND END



5104 EXTRACTION TWEEZER ALBRECHT 16cm POINTED END

FORCEPS | EXTRACTION/RIDGE COMPRESSING FORCEPS


4972S X EXTRACTION FORCEPS SPECIAL MOLAR SET



4973
22 L-R/U-L
SPECIAL MOLAR



4972
22 L-L/U-R
SPECIAL MOLAR



4994PS X ROOT EXTRACTION FORCEPS SET



4997T 51SR
SPECIAL ROOT
UPPER



4994T 47SR
SPECIAL ROOT
LOWER

5033 DR JAMES RIDGE COMPRESSING FORCEP



**5040
ATRAUMATIC UNIVERSAL EXTRACTION FORCEP SET**



4974P
F1P
UPPER
UNIVERSAL



4975P
F2P
LOWER
UNIVERSAL



4976P
F3P
UPPER
ANTERIOR

FORCEPS | EXTRACTION FORCEPS



4977P
F4P
LOWER
ANTERIOR



4979P
F6P
UPPER
MOLAR



4980CP
F7CP
LOWER
MOLAR

FORCEPS | EXTRACTION FORCEPS



FORCEPS | EXTRACTION FORCEPS

5042 FULCRUM FORCEPS KIT



5042-1 FULCRUM EXTRACTION FORCEP LEFT



5042-2 FULCRUM EXTRACTION FORCEP RIGHT

5042-P1
FULCRUM BLUE
THICK PAD 12PCS



5042-P2
FULCRUM WHITE
THIN PAD 12PCS



4033M-1 POWER FLEX TWIST PERIOTOME/ROOT PICK ANTERIOR LEFT S/E



4033M-2 POWER FLEX TWIST PERIOTOME/ROOT PICK ANTERIOR RIGHT S/E



4033M-3 POWER FLEX TWIST PERIOTOME/ROOT PICK POSTERIOR D/E



OPERATIVE



G

CHISELS


Chisels are used for dislodging large bridges of calculus in the interproximal areas. They are also used for cavity preparation, forming line angles on anterior areas and producing proper level on enamel margins.



2101 H8/9 (10-7-14)



2102 H13/14 (20-9-14)



2106 H53/54 (10-6-12)



2103 H15/16 (15-8-14)




2107 H14LR (15-8-14)



2104 H17/18 (10-6-14)



2108 H15LR (15-8-14)



2105 H51/52 (15-8-12)



2111 C1/2 (20-15-3)

OPERATIVE | CHISELS



2112 C3/4 (11-5-3)



2113 C5/6 (15-15-3)



2114 C7/10



2116 C8/9



2118 C40/41 (18-10-6)



2119 MT26 (13-95-8-14)



2120 MT27 (13-80-8-14)



2121 MT28 (10-95-7-14)





2122 MT29 (10-80-7-14)




2123 MT77/78 (15-95-8-12)

EXCAVATORS

Excavators are used to prepare cavities for restorative procedures. They can also be used to remove crown pulp during pulp opening, or to scrape away necrotic tissue in cases of gingival papilla necrosis.




OPERATIVE | EXCAVATORS






PLUGGERS/ CONDENSERS

Pluggers and condensers are used to place filling material into cavities and shape it evenly so it doesn't pop out. They're also used to pack gutta percha into root canals during root canal procedures.



OPERATIVE | PLUGGERS/CONDENSERS



CARVERS

Carvers are used in restorative procedures to remove excess material, carve anatomical features, and contour plastic dental fillings. They can also be used to shape dental wax during large fillings.



2501 IPC



2502 IPCOA



2506 IPC-1




2502T IPCOA TITANIUM



2507 DISCOIDS VIGNON




2503 CVHL1/2



2508 DISCOIDS CD89/92 4EZ



2504 CVHL3S



2508E DISCOIDS CD89/92 6EZ



2505 CVHL3



2509 DISCOIDS CD1/2 4EZ

OPERATIVE | CARVERS



2509E DISCOIDS CD1/2 6EZ



2510 DISCOIDS CD3/6



2511 DISCOIDS CD4/5



2512 DISCOIDS CD90



2513 DISCOIDS WACD



2515 TAN3



2523 WARD2



2608 E2/3



2620 VEHE

CARVING KNIVES

Dental carving knives are used to shape and contour dental materials such as amalgam, composite resin or other restorative substances. They are also used to trim excess filling materials.



2702 21



2707 14



2708 7



2711 30



2712 GF1169

BURNISHERS

Burnishers are used to smooth and polish amalgam and composite fillings, and to remove excess material. They are often used at the end of dental procedures or surgeries.




OPERATIVE | BURNISHERS



OPERATIVE | BURNISHERS



OPERATIVE | BURNISHERS



2822 B33L 2.6mm BALL x 14.5mm TITANIUM



2823 B32M 1.9mm BALL x 12mm TITANIUM



2824 B31S 1.9mm BALL x 10mm TITANIUM



2825T RU1 TITANIUM



2826T RU2 TITANIUM



2827T RU3 TITANIUM



2840 ACORN STR SERRATION 3.8/3.2mm 6EZ

PLASTIC FILLING

Plastic filling instruments, also known as placement instruments, are used to place and shape restorative materials into cavities and during other dental procedures. They are used to smooth out rough edges and irregularities to create a more natural-looking restoration.



2901 PFI T1



2901T PFI T1 TITANIUM



2904T PFI T4 TITANIUM




2902 PFI T2



2905 M1



2902T PFI T2 TITANIUM



2906 M3



2903 PFI T3



2907 M4



2903T PFI T3 TITANIUM



2909 W3



2904 PFI T4



2910 TUFTS4 4EZ

OPERATIVE | PLASTIC FILLING



2910E TUFTS4 6EZ



2910T TUFTS4 TITANIUM



2911 IPC



2912 PFI4/5



2913 PFI6 4EZ



2915 1



2916 B1




2917 H1



2919 WDS3



2920 WDS2



2920T WDS2 TITANIUM




2921T PFI T1



2922T PFI T1




2923T PFI T1



2931 SUPER PLUGGER SMALL/MEDIUM



3001 PFI8A 4EZ



3001E PFI8A 6EZ



3001T PFI8A TITANIUM



3002 PFI11



3003 PFI179



3002T PFI11 TITANIUM



3004 PFIG




3006 TRICO1/2



3007 PFIW-1



3011 IPCOA



3012 PFIH2



3016 PF3R



3017 PF2R

GINGIVAL RETRACTORS/ CORD PACKERS

Gingival retractors and cord packers are used to place and pack retraction cords into the sulcus of a tooth. They create space between the tooth and the gum tissue for taking dental measurements and impressions.



3102 GRM1



3103 GRM2



3108 GCP113



3104 GRM3



3108S GCP113 SERRATED



3105 GRM4 6EZ



3109 GCSI1



3106 BN1



3109S GCSI1 SERRATED



3106S BN1 SERRATED



3110 GCP7



3107 GCPS6 NON-SERRATED



3110S GCP7 SERRATED

CORD PACKERS/ PKT WAXING

Dental cord packers and PKT (packable, kind of soft) waxing are tools used primarily in restorative dentistry, particularly during procedures involving impression making and crown preparations. PKT waxing refers to a technique using specific types of waxes for creating patterns or models in restorative procedures.



3121 GCPS171 NON SERRATED



3121N GCPS171



3203 PKT3R



3121NS GCPS171S WITH SERRATION



3204 PKT3



3121S GCPS171S WITH SERRATION



3206 PKT5



3122 PASCAL



3211 FOGEL-B CCCI-B



3122S PASCAL SERRATED



3211T FOGEL-B CCCI-B TITANIUM



3125 GAUZE PACKER



3212 FOGEL-A CCCI-A

OPERATIVE | CORD PACKERS/PKT WAXING



3212T FOGEL-A CCCI-A TITANIUM



3213T FOGEL-C CCCI-C TITANIUM



3214T FOGEL-D CCCI-D TITANIUM



3215T FOGEL-E CCCI-E TITANIUM



3216T FOGEL-F CCCI-F TITANIUM




3217T FOGEL-G CCCI-G TITANIUM



3218T FOGEL-H CCCI-H TITANIUM

AMALGAM CARRIERS

Amalgam carriers are used to transport, collect, and dispense amalgam or composites into cavities during restorative procedures. They allow for controlled placement of the amalgam into the cavity, helping to minimize waste and ensure that the material fills the space accurately.



4308 S/M 2.0/2.8mm



4309 M/L 2.8/3.4mm



4310 M/J 2.8/3.9mm



4310EZ M/J 2.8/3.9mm EZ REPLACE TIP



4311 AMALGAM WELL SMALL

OPERATIVE TITANIUM

Titanium coated operative instruments that come in a variety of sizing. The titanium coat offers a natural non-stick surface so that you can sculpt material with ease.




2404T TAN3 TITANIUM



2407T WARDS3 TITANIUM



2421T M/L TITANIUM



2413T OREGON1 TITANIUM




2422T KTO TITANIUM



2413TR OREGON1 ROUND TITANIUM




2501T IPC TITANIUM



2415T OREGON3 TITANIUM




2502T IPCOA TITANIUM




2415TR OREGON3 ROUND TITANIUM



2801T B18 TITANIUM



2420T S/M TITANIUM



2802T B26/27S TITANIUM


OPERATIVE | OPERATIVE TITANIUM



2803T B27/29 TITANIUM



2903T PFI T3 TITANIUM



2806T B21 TITANIUM




2904T PFI T4 TITANIUM



2807T B21L TITANIUM




2910T TUFTS4 TITANIUM



2820T B28/29 TITANIUM



2920T WDS2 TITANIUM



2825T RU1 TITANIUM



2921T PFI CM-11 TITANIUM



2826T RU2 TITANIUM



2922T PFI M-12 TITANIUM



2827T RU3 TITANIUM



2923T PFI T-5 TITANIUM



2901T PFI T1 TITANIUM



3001T PF18A TITANIUM



2902T PFI T2 TITANIUM



3002T PF111 TITANIUM

MATRIX RETAINERS

Matrix retainers are used in restorative dentistry to hold a matrix band in place around a tooth during the placement of dental materials, such as composite resin or amalgam. This creates a mold or form to contain the restorative material while it sets.



6204 IVORY 8N



6207 SIQVELAND WIDE



6208 TOFFLEMIRE



6208N
TOFFLEMIRE
LARGE GRIP

6208NM
TOFFLEMIRE
LARGE GRIP
MOLAR BANDS
25 PCS



6210 TOFFLEMIRE JUNIOR



6212 IVORY



6213 IVORY



6214 IVORY SEPARATOR

MATRIX FORCEPS

Matrix Forceps are dental instruments that are utilized to place, remove, and maneuver matrix bands and rings.



7236 IVORY LIGHT WEIGHT 17cm



6216 SECTIONAL FORCEPS



5416 CASTROVIEJO MATRIX CVD 13cm



5531N MATRIX PLACEMENT & REMOVAL CVD

TESTIMONIALS



"Power Dental Group consistently delivers cutting-edge technology that enhances our workflow. Their reliable products have become essential in streamlining our operations, making each day more productive."

- Dr. Jesus Botello OMFS, M.D.



"In all my years as a clinician, I've never encountered instruments that make such a noticeable difference in my workflow until I started using Power Dental Group's instruments. They're not just durable; they make my procedures smoother and more efficient. I can move through surgeries with greater ease and confidence, knowing these instruments will perform without fail. It's clear they've been designed with innovation and precision in mind. For anyone serious about improving their practice, Power Dental Group is the only choice I'd recommend."

- Dr. Mark Lin B.Sc., DDS., M.Sc., F.R.C.D.(C), CEO and Founder of T.I.D.E. Canada



"Power Dental Group crafts the finest and most durable dental and surgical instruments in the world. Tested through thousands of rigorous surgeries, our instruments consistently deliver flawless performance. For clinicians who demand unparalleled luxury, innovation, and uncompromising quality, trust Power Dental Group to be an extension of your clinical team."

- Dr. Ady Palti DDS, Chairman of the ICOI



"I used your Crown Removers to remove a zirconia crown on 19 to do an endodontic retreat. It worked so well that it surprised me and my assistant. Thanks for making it. This saves me valuable chair time and hundreds of dollars in lab bills."

- Dr. Shawn Lottier DDS, MAGD, DICOI



"Power Dental Group has continued to push the envelope in innovation with their patented instrumentation. They make my life easier by having high integrity products that radically improve efficiency in my day-to-day."

- Dr. Eric Park DDS



"Power Dental Group has revolutionized our daily tasks with their advanced equipment. Their dedication to high standards ensures we can perform at our best, day in and day out."

- Dr. Mark Kwon DMD



"I've worked with a lot of different instruments over the years, but Power Dental Group's are a game changer. Their instruments are incredibly efficient - it's like they were made with a deep understanding of what a surgeon really needs. They help me move faster, stay precise, and get better outcomes for my patients. Honestly, I can't imagine going back to anything else."

- Dr. David Wong DDS, Periodontist



ENDODONTIC



H

PLUGGERS & CONDENSERS

Pluggers and condensers are used to place filling material into cavities and shape it evenly so it doesn't pop out. They're also used to pack gutta percha into root canals during root canal procedures.



1026 DG 16 3EZ



3305 L1



3306 L2



3307 1/3



3308 RCP5/7



3309 RCP9/11

ENDODONTIC | PLUGGERS & CONDENSERS



3312 GLICK1



3331 GLICK2



3313 RCA8A



3314 RCA8.5A



3315 RCA9A



3316 RCA9.5A



3317 RCA10A



3318 RCA10.5A

ENDODONTIC | PLUGGERS & CONDENSERS



3319 RCA11A



3320 RCA11.5A



3321 RCA12A



3322 RCP8P



3323 RCP8.5P



3324 RCP9P



3325 RCP9.5P



3326 RCP10P



3327 RCP10.5P

ENDODONTIC | PLUGGERS & CONDENSERS



3328 RCP11P



3329 RCP11.5P



3330 RCP12P



3413 GP2



3414 HCO



3415 HC00



3416 MA57



3417 GP3



3418 DG16

3332D

NiTi ENDO PLUGGER SET OF 3



3332-1 NiTi ENDO PLUGGER 0.30/0.75



3332-2 NiTi ENDO PLUGGER 0.40/0.80



3332-3 NiTi ENDO PLUGGER 0.75/1.30

SPREADERS & HEAT CARRIERS

Spreaders are used to pack root canal filling material into the sides of a canal. Heat carriers are a type of dental spreader with enhanced heat retention properties for gutta-percha compaction.



ENDODONTIC | SPREADERS & HEAT CARRIERS



3405 GP1



3406 W1S



3408 D11



3410 D11T



3411 D11TS

RUBBER DAM CLAMP ACCESSORIES

Rubber dam clamp accessories are tools and components used in conjunction with rubber dams during dental procedures. Rubber dams are used to isolate the treatment area, improve visibility, and maintain a dry field.



7237 BREWER 17.5cm



7239 AINSWORTH 17cm



7238 STOKES 17.5cm



7240 IVORY 17cm



7238N RING PLACEMENT 15cm



7241S RUBBER DAM FRAME 9cm
7241M RUBBER DAM FRAME 13cm



CROWN INSTRUMENTS



CROWN INSTRUMENTS

Crown instruments are specialized tools used in the preparation, placement, and adjustment of dental crowns. These instruments help ensure that crowns fit well, look natural, and function properly. By ensuring precision and efficiency throughout the process, these instruments contribute significantly to the overall quality and efficiency of dental restorations.



4602
CROWN
SCISSOR CVD
10.5cm
SERRATED



4602T
CROWN
SCISSOR
CVD 10.5cm
SERRATED T/C



4603
CROWN
SCISSOR
UNIVERSAL
12cm



4615
CROWN
SCISSOR
STR

CROWN INSTRUMENTS | CROWN INSTRUMENTS




4615T
CROWN
SCISSOR
STR T/C



4630
PEDODONTIC
CROWN
SCISSOR
CVD 9cm



4605 CROWN SPREADER CROSS TYPE




4605M CROWN SPREADER FLAT TYPE



4605N CROWN SPREADER
FLAT TYPE SMALL HANDLE



4607T
REPLACEMENT
TIPS 10EA



4607
CROWN
REMOVING PLIERS
W/ SILICONE TIPS



4608
CROWN
SPREADER
FORCEPS



4610
CROWN
GRIPPER T/C



4632
PEDO CROWN
GRIPPER



4706 ARTICULATING PAPER
FORCEPS 16cm STR



4707 MILLER FORCEPS 15cm CVD

CROWN REMOVERS

Crown removers are specialized instruments designed to facilitate the safe and effective removal of dental crowns, whether they are temporary or permanent. They are used to detach crowns from the underlying tooth structure. This is often necessary if the crown is loose, damaged, or requires replacement.

4638 X POWER CROWN REMOVER KIT



4638-R
POWER CROWN
REMOVER - R



4638-L
POWER CROWN
REMOVER - L




4636
SPECIAL
HIGH SPEED
DRILL BUR
2.0mm GOLD



4637
SPECIAL
HIGH SPEED
DRILL BUR
2.0mm SILVER

CROWN INSTRUMENTS | CROWN REMOVERS



CROWN INSTRUMENTS | CROWN REMOVERS



4701 CROWN REMOVER LONG BUMPER




4704 CROWN REMOVER AUTOMATIC




4710T CROWN REMOVER TEETH D/E

IMPRESSION TRAYS

Impression trays are used to hold impression materials (such as alginate, polyvinyl siloxane, or polyether) to create detailed molds of the teeth and gums. They are available in multiple sizes to accommodate different mouth sizes and shapes, ensuring a proper fit for each patient.

		U-1(XS)	U-2(S)	U-3(M)	U-4(L)	U-5(XL)
	A	52mm	54mm	58mm	59mm	63mm
	B	65mm	67mm	72mm	75mm	84mm

		U-1(XS)	U-2(S)	U-3(M)	U-4(L)	U-5(XL)
	A	50mm	53mm	57mm	60mm	66mm
	B	65mm	68mm	75mm	78mm	84mm



7301UXS UPPER XS
7301US UPPER S
7301UM UPPER M
7301UL UPPER L
7301UXL UPPER XL



7301LXS LOWER XS
7301LS LOWER S
7301LM LOWER M
7301LL LOWER L
7301LXL LOWER XL





7302UXS UPPER XS
7302US UPPER S
7302UM UPPER M
7302UL UPPER L
7302UXL UPPER XL



7302LXS LOWER XS
7302LS LOWER S
7302LM LOWER M
7302LL LOWER L
7302LXL LOWER XL

CROWN INSTRUMENTS | IMPRESSION TRAYS

	UPPER					
		U-1(XS)	U-2(S)	U-3(M)	U-4(L)	U-5(XL)
	A	56mm	58mm	63mm	64mm	68mm
	B	62mm	64mm	70mm	77mm	78mm

		U-1(XS)	U-2(S)	U-3(M)	U-4(L)	U-5(XL)
	A	53mm	56mm	60mm	62mm	63mm
	B	62mm	65mm	72mm	75mm	83mm



7303UXS UPPER XS
7303US UPPER S
7303UM UPPER M
7303UL UPPER L
7303UXL UPPER XL



7303LXS LOWER XS
7303LS LOWER S
7303LM LOWER M
7303LL LOWER L
7303LXL LOWER XL



7304UXS UPPER XS
7304US UPPER S
7304UM UPPER M
7304UL UPPER L
7304UXL UPPER XL



7304LS LOWER S
7304LM LOWER M
7304LL LOWER L
7304LXL LOWER XL



7308U
UPPER
65X39mm



7308L
LOWER
64X39mm



7309L
UPPER LEFT
LOWER RIGHT
60mm



7309R
UPPER RIGHT
LOWER LEFT
60mm

7310
ADJUSTABLE
PERFORATED



7312
BITE
CHECA
TRAY

7314

DISPOSABLE FISH BONE TRAY SET OF 5



7314-1 POSTERIOR QUADRANT 15mm WALL 35EA

7314-2 POSTERIOR SIDE 14mm WALL 35EA

7314-3 POSTERIOR HALF SIDE 7mm WALL 50EA

7314-4 POSTERIOR SIDELESS 4mm WALL 50EA

7314-5 ANTERIOR 20mm WALL 50EA



7412

ENDODONTIC BOX
10cmX4.5cmX4.2cm



7421

SPIRIT LAMP



7425

BONE TRAY SQUARE
65x58x13mm



7426

BONE TRAY BOWL
58mm DIAx28mm DEPTH



ORTHODONTIC



J

CUTTERS & PLIERS

Cutters and pliers are used in various dental procedures for cutting, shaping, and manipulating dental materials and tissues. They are used with dental handpieces to cut through tooth structure or dental materials during preparations for restorations.



5601
DISTAL END
CUTTER T/C
SAFTY HOLD
12.5cm



5601F
DISTAL END
FLUSH
CUTTER
T/C SAFTY
HOLD 13.5cm



5601NXL
DISTAL END
CUTTER T/C
SAFTY HOLD
LONG HANDLE
15cm



5601S
DISTAL END
CUTTER T/C
SAFTY HOLD
SLIM
12.5cm



5602
DISTAL END
CUTTER T/C
SAFTY HOLD
11.5cm



5602HN
DISTAL END
FLUSH
CUTTER T/C
SAFTY HOLD
12cm



5603N2
DISTAL END
CUTTER T/C
SLIM 15°
ANGLED 13cm
J-1



5604
WIRE
CUTTER
STR T/C
12cm



5607K
WIRE CUTTER
STR THIN
SOFT T/C
12.5cm



5609
WIRE CUTTER
STR THIN
SOFT T/C 15°
ANGLED 11.5cm



5610
PLIER
WEINGART
THICKER



5611K
PLIER
WEINGART
T/C SLIM

ORTHODONTIC | CUTTERS & PLIERS



5611N
PLIER
WEINGART
SLIM O



5612
PLIER
HOW
STR



5617
PLIER
ADHESIVE
REMOVER



5617MT
REPLACEMENT
METAL TIPS



5617T
REPLACEMENT
PLASTIC TIPS 3EA



5618
PLIER
BAND
REMOVER



5618T
REPLACEMENT
PLASTIC TIPS 3EA



5619
PLIER
BAND
REMOVER



5620
PLIER
BRACKET
REMOVER
STR



5621
PLIER
CONTOURING



5622
PLIER
HOW
THIN



5624
WIRE CUTTER
STR HARD T/C



5632
WIRE CUTTER
SOFT T/C

PLIERS

Pliers are used to manipulate wires, place brackets, and adjust appliances with precision. The jaws can be straight, curved, or angled, depending on the specific type of plier.



5701
WALDSACHE



5703
NANCE



5704
ADERER



5705
GOSLEE



5706
MODELL
MARBURG
SERRATED



5709
ARROW
CLASP
BENDING



5710
ARROW
CLASP
FORMING



5712
PEESO



5712S
PEESO
SHORTER
BEAK



5714
LANGENBEC



5716
JARABEK



5718
BIRD BEAK W/
SOFT WIRE
CUTTER 12cm

ORTHODONTIC | PLIERS



5718J2
BIRD BEAK
W/ SOFT
WIRE
CUTTER 2
GROOVE



5718L
BIRD LONG
BEAK NO
GROOVE
13.5cm



5718LN
BIRD LONG
BEAK
1 GROOVE
13cm



5718US
BIRD
BEAK
12.5cm



5721
SCHWARZ
NEW



5723
3-PRONG
REGULAR
12.5cm



5723S
3-PRONG
SMALL
12cm



5724
TWEED
12cm



5725
LOOP FORMING/
TWEED PLIER
O'BRIEN 12.5cm



5725J
LOOP FORMING/
TWEED PLIER
3 BARREL W/ SOFT
WIRE CUTTER



5726
WIRE
BANDING



5727
YOUNG



5728
CLASP
FORMING



5729
ADAMS



5730
FREVERT



5731
NANCE
LOOP
BENDING



5732
HOLLOW
CHOP ARCH
FORMING 3
GROOVES



5733
STEP
BANDING
0.75mm

ORTHODONTIC | PLIERS



5734
STEP
BANDING
1mm



5736
PUSHING
& TIGHT
T/C 12cm



5739
STEP
BANDING
0.25mm



5740
STEP
BANDING
0.5mm



5743
WIRE
BENDING
T/C



5744
CONTOURING
HOLLOW
CHOP



5745
CINCH
BACK



5748
CROWN
& BAND
CRIMPING



5749
BAND
CRIMPING



5750
BAND
REONALD



5754
TWEED



5765
LIGHT
WIRE
BENDING



5770
PEDO
CRIMPING



5771
BENDING



5772
OPTICAL
BENDING



5774
CRIMPABLE
HOOK



5776
TORQUE
BANDING
PLIER

ORTHODONTIC | PLIERS



5777
BRACKETS
REMOVING
FOR CERAMIC



5778
CRIMPABLE
HOOK
PLACEMENT




5780
V-SHAPE
BENDING



5781
LINGUAL
ARCH
FORMING



5782
UTILITY



5741 PLIER RACK 22cm



ORTHODONTIC INSTRUMENTS

Orthodontic instruments are used to diagnose, treat, and manage dental and facial irregularities. These instruments play a crucial role in various procedures, from initial assessments to the final adjustments of orthodontic appliances.



5746 HEIGHT GAUGE 0.022

5751 HEIGHT GAUGE 0.018



5817
BRACKET
POSITIONING
GAUGE



5755 DISC GUARD SMALL 22mm
OUT 20mm IN $\pm 3\%$

5756 DISC GUARD LARGE 25mm
OUT 23mm IN $\pm 3\%$



J-7 **5775** TORQUE BENDING PLIER TWIN HEAD WITH KEY



5753
CONVERTIBLE
CAP
REMOVER



5809
SEPARATING
15.5cm



5809S
SEPARATING
15.5cm
S CURVE



5810
MATHIEU
PLIERS
14cm



5810J
LIGATURE
PLACING
PLIER HOOK
SHAPE



5812
MATHIEU
PLIERS
14cm
HOLLOW BACK



5812C
MATHIEU
PLIERS
14cm
HOLLOW
BACK CVD



5814
MOSQUITO
PLIER
STR



5818
LIGATURE
TYING
PLIERS



5842
WIRE
TWISTER
FORCEP



5807 BRACKET HOLDER 13.5cm



5808 BRACKET HOLDER 14cm



5816 BRACKET PLACING 12.5cm



5816L BRACKET PLACING DIRECT BOND 203-O



5816S BRACKET PLACING KEAT MICRO 12.5cm



5816T BRACKET PLACING DIRECT BOND WITH BRACKET ALIGNER



5819 BRACKET HOLDING 15cm



5843 CHAIN HOLDER 1



5844 CHAIN HOLDER 2



5845 CHAIN HOLDER 3



5862 LOTUS TWEEZER SLB 202-O



5863 KEAT MICRO 12.5cm 209-O



5864 LOTUS PLUS ARCHWIRE 211-O



5865 LOTUS OPEN/CLOSING 210-O



5866 PROSELECT/SLOT ALIGNER 200-322DF



5867 DELUX BRACKET HOLDER

ORTHODONTIC | ORTHODONTIC INSTRUMENTS



5801 BAND SEAT & REMOVER



5801N BAND SEAT & REMOVER NEW



5804 ELASTIC INSTRUMENT 15.5cm



5833 BAND SEATER/SCALER



5840 U15/FILE



5841 POINT/FILE



5850 LIGATURE DIRECTOR 200-201



5851 LIGATURE DIRECTOR 201-501



5852 LIGATURE DIRECTOR 201-502 SCALER/BAND THINNER

ORTHODONTIC | ORTHODONTIC INSTRUMENTS



5853 SHURE BAND SEATER 200-208 SCALER/BAND THICKER



5854 PICKER/EXPLORER



5855 BAND PUSHER/SCALER 2100-O



5856 DIRECT BAND ADHESIVE REMOVER 2102-O



5857 LIGATURE DIRECTOR 2104-O



5858 LIGATURE DIRECTOR/EXPLORER 2106-O



5859 DISTAL BANDER 0.018" AND 0.030"



5868 PROSELECT BRACKETIZER



5869 BAND PUSHER/SCALER 201-500



5869T BAND PUSHER/SCALER T



5870 LIGATURE REMOVER/DIRECTOR 201-503



5871 ELASTIC REMOVER 200-109



5872 BAND PUSHER/SCALER



5873 ELASTIC REMOVER/CP15



5874 EXP23/PLG2.7 SERRATION



5875 EXP23/LIGATURE



5876 EXP23/LIGATURE



5877 BAND SEAT OBP S/E



5877D BAND SEAT OBP D/E



5878 LIGATURE OLD1/U15



5879 LIGATURE OLD1 S/E



5880 LIGATURE OLD1/BAND SEAT OBP



5881 J1/LIGATURE OLD1



5882 U15/PLG 2.8 SERRATION



5883 LIGATURE OLD PLG 0.7 SERRATION



5884 LIGATURE OLD PLG 1.5 SERRATION



5885 LIGATURE OLD PLG 1.7 SERRATION

ORTHODONTIC | ORTHODONTIC INSTRUMENTS



5886 LIGATURE OLD PLG 2.5 SERRATION



5887 PLG 1.1/2.8 SERRATION



5889 EXP23 1/8 / EXP6 5.3mm 1/8



5890 EXP23/BAND PUSHER



5893 PICKER/LIGATURE DIRECTOR



5894 BAND SEAT OBP/H5



5895 BAND SEAT OBP/U15



5896 BAND PUSHER/LIGATURE



5769 DONTRIX GAUGE
ELASTOMETRICS TENSION FORCE STRESS 16oz (0-400g)



5773 BENDING DRIVER



5805 BAND PUSHER 14cm



5806 BAND SETTER 15cm



5820 BAND SEATING 13cm



5832 SOLDERING TWEEZER



5860 NITI ICE 2000-O



5861 POSITIVE POSITIONER 2222-O



5891 SHORTY TWIST TYING LIGATURE



GENERAL INSTRUMENTS



K

WAX & MODELING CARVERS

Wax and modeling carvers are used to shape and refine wax models, such as those used for dentures, crowns, and bridges, ensuring anatomical accuracy. They can be used to smooth edges and surfaces, enhancing the final aesthetic of dental restorations.



6610 LE CRON 6EZ



6612 BEALE



6701 FAHNENSTOCK 12.5cm

GENERAL INSTRUMENTS | WAX & MODELING CARVERS



6702 LESSMANN 12.5cm



6703 CON SCODELLINO 12.5cm



6704 FAHNENSTOCK 17cm



6705 LESSMANN 17cm



6706 CON SCODELLINO 17cm

SPATULAS

Spatulas are used for mixing dental materials such as cements, impression materials, and composites. They are also used to apply materials to the tooth surface or within a cavity.



6803 DOUBLE END SMALL



6804 DOUBLE END MEDIUM



6805 DOUBLE END LARGE



6806 SINGLE END SMALL



6807 SINGLE END MEDIUM



6808 SINGLE END LARGE

GENERAL INSTRUMENTS | SPATULAS



6812 WOOD HANDLE 21.5cm



6815 WOOD HANDLE 20cm



6818 WOOD HANDLE 21.5cm



6820 LAB KNIFE 14cm



6821 PERIO SPATULA 6mm



6822 PERIO SPATULA 9mm




6830 HOT PLATE SPATULA

SUCTION TUBES

Suction tubes are designed to remove fluids, debris, and other materials from the oral cavity during various dental procedures. They eliminate saliva from the treatment area, keeping the field dry and improving visibility.



GENERAL INSTRUMENTS | SUCTION TUBES



6710E SURGICAL SUCTION 6mm BLACK EDITION



6731 1.5mm

6732 2.5mm

6733 3mm

6734 4mm



6738 2.4mm

6739 3.5mm

6740 4.5mm



6741 2.4mm



6743 4.75mm



6742 3.0mm

TOWEL CLAMPS

Towel clamps are used to hold drapes and towels in place, ensuring a clean and organized working area during dental procedures. They feature a locking mechanism that securely holds towels in place.



7101
BACKHAUS
8cm



7102
BACKHAUS
10.5cm



7105
BACKHAUS
13.5cm

GENERAL INSTRUMENTS | TOWEL CLAMPS



7110
20cm



7112
MAGILL CATHETER
FORCEP
20cm




7113
MAGILL CATHETER
FORCEP
25cm

CASSETTES

Cassettes provide a structured way to organize dental instruments, ensuring that everything is easily accessible and in one place during procedures. By securely holding instruments in place, cassettes help prevent damage during transportation, reducing wear and tear.



7501 CASSETTE 5 HOLD
60x180x17.5mm




7502D CASSETTE DOUBLE HINGE
16 HOLD DOUBLE RACK
114x203x31mm




7502N POWER CASSETTE MEDIUM
148x210x35mm



7502P X CASSETTE OPEN TYPE
205x145x35mm




7502S CASSETTE DOUBLE HINGE
8 HOLD SINGLE TALL RACK
114x203x31mm



7502U CASSETTE DOUBLE HINGE
UTILITY
114x203x31mm

GENERAL INSTRUMENTS | CASSETTES



7503D CASSETTE DOUBLE HINGE
12 HOLD DOUBLE RACK
90x203x31mm



7503N POWER CASSETTE SMALL
95x210x35mm



7503P X CASSETTE OPEN TYPE
205x93x35mm



7503PB X CASSETTE OPEN TYPE
205x93x35mm



7503U CASSETTE
DOUBLE HINGE UTILITY
90x203x31mm




7505-1 CASSETTE CLIPS SCISSOR



7505-2 CASSETTE CLIPS SUCTION

GENERAL INSTRUMENTS | CASSETTES



7511 CASSETTE
SURGICAL/UTILITY
174x203x31mm



7512 X CASSETTE
SURGICAL/UTILITY
275x205x35mm



7530 CASSETTE PIEZO TIPS BUR
HOLDER 6 HOLD



7533 CASSETTE
UTILITY PARTS BOX
80x30x25mm



7534 X CASSETTE
UTILITY 5 HOLD



KITS & SETS



KITS & SETS

KITS & SETS

Crafted from years of clinical understanding and thousands of clinicians' feedback, our kits and sets will provide you with the most effective solutions that will allow you to receive the best results without having to waste your time.

7801

POWER AQUA SINUS LIFT KIT



7802

POWER CRESTAL APPROACH KIT



KITS & SETS

7803

POWER LATERAL APPROACH KIT



7804

BONE SPREADER 2.5mm




7805

BONE TRIMMER KIT




7806

POWER BUR KIT



7808

POWER BONE SPREADING KIT



KITS & SETS

7809

POWER BONE EXPAND & SPREADING KIT



7810

POWER BONE COLLECTING KIT



7811

POWER WATER RINSING SYSTEM



KITS & SETS

7812

POWER DRILL STOP SYSTEM



7813

POWER PROSTHETIC UNIVERSAL DRIVER KIT



7814

POWER EASY BONE EXPANDER KIT



7815

POWER IMPLANT REMOVER KIT



KITS & SETS

7818-1

BONE TACK TITANIUM 3mm 5EA



7818-2

BONE TACK TITANIUM 5mm 5EA



7818-3

BONE TACK CASE 3mm



7818-4

BONE TACK CASE 5mm



7818-7

BONE TACK HOLDER & REMOVER 1



7818-8

BONE TACK HOLDER & REMOVER 2

7818-10

POWER TAC SCREWS 3mm
10EA WITH CASE



7818-11

POWER TAC SCREWS 5mm
10EA WITH CASE



KITS & SETS

7818N

POWER BONE TAC KIT WITH CASE



7819

DRILL GUIDE & STOPPER SET



7822

POWER SCREW & IMPLANT FIXTURE REMOVER KIT



KITS & SETS

7824N

RGF CENTRIFUGE MACHINE



7824-2

PRF TUBES - RED
100/BOX



7824-3

PRF TUBES - YELLOW
100/BOX



7824-4

BLOOD COLLECTION
BUTTERFLY NEEDLE SET 25EA



7824-5

TOURNIQUET
25EA



7824-7

TUBE STAND



KITS & SETS

7824-1 PRF/RGF INSTRUMENTS KIT (10 PCS)



2034S BONE GRAFT PACKER 3/5mm SERRATED



2035T BONE CARRIER TITANIUM



6505 SEMKEN-TAYLORS STR TEETH 1X2 12.5cm



6834 PRF/RGF CARRIER



7824-6 PRF/RGF CASSETTE




5907S
DEAN SUPER CUT 17.5cm



2094
BONE MIXING BOWL



7425
BONE TRAY SQUARE
65x58x13mm



7502U CASSETTE DOUBLE HINGE
UTILITY
114x203x31mm

KITS & SETS

7830

X EASY BONE EXPANDER/SINUS LIFT/ROOT EXTRACTION
WITH UNIVERSAL LATCHET DRIVER SET OF 7



7830A

DR. ADY PALTI X BLACK EDITION IMPLANT SURGERY KIT SET OF 15



KITS & SETS

7854

MAURI TUNNELING KIT 8 INSTRUMENTS WITH CASSETTE



1812M KOR1/2 MICRO BLADE



1953 TUNNELING PROCEDURE BLUE TITANIUM



1954 TUNNELING PROCEDURE LONG BLUE TITANIUM



2044T PPAELA TITANIUM




3121 GCPS171 NON SERRATED




3516T PERIOSTEALS BUSER



5873 ELASTIC REMOVER/CP15



6934 BLACK EDITION FOR MICRO BLADE



7503U CASSETTE
DOUBLE HINGE UTILITY
90x203x31mm


KITS & SETS

7855

MAURI PERIO KIT 18 INSTRUMENTS WITH CASSETTE



1812M KOR1/2 MICRO BLADE



1954 TUNNELING PROCEDURE
LONG BLUE TITANIUM



2038T BONE GRAFT/CARRIER
2-4-6-8-10mm TITANIUM



2044T PPAELA TITANIUM



3121 GCPS171 NON SERRATED



3506T PERIOSTEALS MOLT9



3509T PERIOSTEALS PR3



3516T PERIOSTEALS BUSER



6934 BLACK EDITION FOR MICRO BLADE



6920 6EZ



5873 ELASTIC REMOVER/CP15



6505 SEMKEN-TAYLORS STR TEETH 1X2 12.5cm



6511 SUTURE CORN TISSUE FORCEPS CVD 16cm



5514B
CASTRO
STR T/C 18cm BLACK EDITION



5929S
CASTRO MICRO STR 14cm



6420
MINNESOTA RETRACTOR



6525 GERALD STR 18cm



7512 X CASSETTE
SURGICAL/UTILITY
275x205x35mm

KITS & SETS

7856

MAURI EXTRACTION KIT 5 INSTRUMENTS




3747 BERNARD 3 LEFT



3748 BERNARD 2 STR



3749IN BERNARD IN



3749OUT BERNARD OUT



3750 BERNARD 4 RIGHT

KITS & SETS

9001-ABND ACCU LUX-TOME SPADE SET



3771 ACCU LUX-TOME SPADE



3790 ACCU BERNARD SPADE OFFSET IN



3791 ACCU BERNARD SPADE OFFSET OUT



3792 ACCU BERNARD SPADE LEFT



3793 ACCU BERNARD SPADE RIGHT



3794 ACCU LUX-TOME T-BAR



7503U CASSETTE
DOUBLE HINGE UTILITY
90x203x31mm

KITS & SETS

9002-ADSI ADVANCE SINUS SET



1970 SINUS BONE GRAFT 3.87SP/2.8mm (2-4-5-6-8-10mm)



1971 SINUS COBRA 3.5mm L/R



1972 SINUS COBRA IN/OUT



2069 MUSHROOM 2.8mm/3.8mm (2-4-5-6-8-10)



2070 MUSHROOM 4.8mm/5.8mm (2-4-5-6-8-10)



7503U CASSETTE
DOUBLE HINGE UTILITY
90x203x31mm

KITS & SETS

9003-AMOP ALN MICRO SURGERY KIT 26PCS



1103E CP11 (3-6-8-11mm)



1505 YC7/8



1812 KOR1/2



1812M KOR1/2 MICRO BLADE



1901 PR1/2



2011 OCH2



2043T PPAEL TITANIUM



2044T PPAELA TITANIUM



3509 PR3




3515 BENQUE



3516 BUSER



5302C MICRO-HALSTED
12.5cm CVD



5505TN MICRO CRILLE-WOOD
STR T/C 15cm



5505B CRILE-WOOD
BLACK EDITION STR 15cm



5513B CASTROVIEJO
BLACK EDITION STR 16cm



5516B CASTROVIEJO
BLACK EDITION CVD 18cm



5911S GOLDMAN-FOX
SUPER CUT STR 12.5cm



6511 SUTURE CORN TISSUE
FORCEPS CVD 16cm



6511N SUTURE CORN TISSUE FORCEPS
MICRO SUTURE CVD 16cm



6525C GERALD MICRO 18cm




6733 SUCTION TUBE 3mm



6905 5A ROUND STR



6911 MIRROR HANDLE 6EZ C/S



6951D MIRROR #5 24.8mm C/S



7018N SYRINGE 1.8ml



7512 X CASSETTE
SURGICAL/UTILITY
275x205x35mm

KITS & SETS

9004-ATRA ATRAUMATIC EXTRACTION KIT

1105 CP15UNC

1923T PT2 TITANIUM

1924T PT1 TITANIUM

2043T PPAEL TITANIUM

3506T PERIOSTEALS MOLT9

3752 3mm STR

3755 3mm CVD

4107 LUCAS87

4222 BONE FILE 21

4411 BONE RONGEURS
MICRO FRIEDMAN 15 14.5cm

4974 F1 UPPER UNIVERSAL

4975 F2 LOWER UNIVERSAL

4994 46XL LOWER ROOT

5505B CRILE-WOOD
BLACK EDITION STR 15cm

6010

IRIS CVD
T/C 11.5cm

6501 ADSON STR 12cm

6501A ADSON CVD TEETH 12cm

6517 LONDON-COLLEGE

6715 2.5mm

6905 5A ROUND STR

6911 MIRROR HANDLE 6EZ C/S

6951D MIRROR #5 24.8mm C/S

7018N SYRINGE 1.8ml

7505-1 CASSETTE
CLIPS SCISSOR

7505-2 CASSETTE
CLIPS SUCTION

7512 X CASSETTE
SURGICAL/UTILITY
275x205x35mm

KITS & SETS

9005-BACL ACCU-LUXATING SET 8PCS



3769 ACCU LUX-TOME STR SERRATED TITANIUM



3770 ACCU LUX-TOME CVD SERRATED TITANIUM



3771 ACCU LUX-TOME SPADE SERRATED TITANIUM



3772 ACCU LUX-TOME DISTAL SERRATED TITANIUM



3773 ACCU LUX-TOME MESIAL SERRATED TITANIUM



3774 ACCU LUX-TOME LARGE DISTAL SERRATED TITANIUM



3775 ACCU LUX-TOME LARGE MESIAL SERRATED TITANIUM



7503U CASSETTE
DOUBLE HINGE UTILITY
90x203x31mm

KITS & SETS

9006-BAFO

ATRAUMATIC EXTRACTION FORCEP SET 6EA



4974
F1 UPPER
UNIVERSAL



4975
F2 LOWER
UNIVERSAL



4976
F3 UPPER
ANTERIOR



4977
F4 MODIFIED
LOWER ANTERIOR



4979
F6 UPPER
MOLARS



4980C
F7C LOWER
3RD

KITS & SETS

9007-BGRA BONE GRAFTING KIT



1816D KM3/4



2033NT BONE GRAFT PACKER 3.3/4mm TITANIUM



2038T BONE GRAFT/CARRIER 2-4-6-8-10mm TITANIUM



3513T ALLEN9 TITANIUM




6502 ADSON STR TEETH 1x2 12cm

5505B
CRILE-WOOD
BLACK EDITION
STR 15cm




5917B
LA GRANGE
SUPER CUT
BLACK EDITION
CVD 11.5cm



7503D CASSETTE DOUBLE HINGE
12 HOLD DOUBLE RACK
90x203x31mm



6923
SCALPEL HANDLE
180° ADJUSTABLE
16.2cm



KITS & SETS

9008-BIBT PERIOSTEAL ELEVATORS SET



2001T 13K/TG 3.5/3mm TITANIUM



2002T 3K/LR 4.5/3.5mm TITANIUM



3506T PERIOSTEALS MOLT9



3509T PERIOSTEALS PR3



3513T ALLEN9 TITANIUM



3514T 24G TITANIUM



3516T BUSER TITANIUM




3521T PERIOSTEALS HF20



3526T PERIOSTEALS WOODSON-1 MODIFY



3534T PERIOSTEALS PERIO SURGICAL



7502S CASSETTE DOUBLE HINGE
8 HOLD SINGLE TALL RACK
114x203x31mm

KITS & SETS

9009-BOSTC OSTEOTOME SET 8PCS



2050 OSTEOTOME 3.2mm CVD (8-10-13-15-18mm)



2051 OSTEOTOME 4.0mm CVD (8-10-13-15-18mm)



2052 OSTEOTOME 1.8mm CVD (8-10-13-15-18mm)



2053 OSTEOTOME 2.5mm CVD (8-10-13-15-18mm)




2054 OSTEOTOME 3.0mm CVD (8-10-13-15-18mm)



2055 OSTEOTOME 4.0mm CVD (8-10-13-15-18mm)



2056 OSTEOTOME R0.38mm CVD (8-10-13-15-18mm)



7502U CASSETTE DOUBLE HINGE
UTILITY
114x203x31mm

KITS & SETS

9010-BOSTM


OSTEOTOME LONG OFFSET CONVEX SET 6PCS



2095R OSTEOTOME 2.0mm (4-6-8-10-13-16-18-20-23-26mm)
LONG OFFSET CONVEX



2096R OSTEOTOME 3.5mm (4-6-8-10-13-16-18-20-23-26mm)
LONG OFFSET CONVEX



2097R OSTEOTOME 4.3mm (4-6-8-10-13-16-18-20-23-26mm)
LONG OFFSET CONVEX



2098R OSTEOTOME 5.0mm (4-6-8-10-13-16-18-20-23-26mm)
LONG OFFSET CONVEX



2099R OSTEOTOME 5.5mm (4-6-8-10-13-16-18-20-23-26mm)
LONG OFFSET CONVEX



7503S CASSETTE DOUBLE HINGE
6 HOLD SINGLE TALL RACK
90x203x31mm

KITS & SETS

9011-BOSTN

OSTEOTOME LONG OFFSET CONCAVE SET 6PCS



2095 OSTEOTOME 2.0mm (4-6-8-10-13-16-18-20-23-26mm)
LONG OFFSET CONCAVE



2096 OSTEOTOME 3.5mm (4-6-8-10-13-16-18-20-23-26mm)
LONG OFFSET CONCAVE



2097 OSTEOTOME 4.3mm (4-6-8-10-13-16-18-20-23-26mm)
LONG OFFSET CONCAVE



2098 OSTEOTOME 5.0mm (4-6-8-10-13-16-18-20-23-26mm)
LONG OFFSET CONCAVE



2099 OSTEOTOME 5.5mm (4-6-8-10-13-16-18-20-23-26mm)
LONG OFFSET CONCAVE



7503S CASSETTE DOUBLE HINGE
6 HOLD SINGLE TALL RACK
90x203x31mm

KITS & SETS

9012-BOSTP

OSTEOTOME LONG STRAIGHT CONVEX SET 6PCS



2095SR OSTEOTOME 2.0mm (4-6-8-10-13-16-18-20-23-26mm)
LONG CONVEX



2096SR OSTEOTOME 3.5mm (4-6-8-10-13-16-18-20-23-26mm)
LONG CONVEX



2097SR OSTEOTOME 4.3mm (4-6-8-10-13-16-18-20-23-26mm)
LONG CONVEX



2098SR OSTEOTOME 5.0mm (4-6-8-10-13-16-18-20-23-26mm)
LONG CONVEX



2099SR OSTEOTOME 5.5mm (4-6-8-10-13-16-18-20-23-26mm)
LONG CONVEX



7503S CASSETTE DOUBLE HINGE
6 HOLD SINGLE TALL RACK
90x203x31mm

KITS & SETS

9013-BOSTS OSTEOTOME LONG SET 8PCS



2050S OSTEOTOME 3.2mm STR (8-10-13-15-18mm)



2051S OSTEOTOME 4.0mm STR (8-10-13-15-18mm)



2052S OSTEOTOME 1.8mm STR (8-10-13-15-18mm)



2053S OSTEOTOME 2.5mm STR (8-10-13-15-18mm)




2054S OSTEOTOME 3.0mm STR (8-10-13-15-18mm)



2055S OSTEOTOME 4.0mm STR (8-10-13-15-18mm)



2056S OSTEOTOME R0.38mm STR (8-10-13-15-18mm)



7502U CASSETTE DOUBLE HINGE
UTILITY
114x203x31mm

KITS & SETS

9014-BPT PERIOTOME SET 5PCS



KITS & SETS

9015-BSINL SINUS LIFT SET 10EA



2021 M6577



2023 SINUS LIFT 1



2024 SINUS LIFT 2



2025 SINUS LIFT 3



2026 SINUS LIFT 4



2027 SINUS LIFT 5



2062 SINUS LIFT 6



2063 SINUS LIFT 7



2069 MUSHROOM 2.8mm/3.8mm (2-4-5-6-8-10)



2070 MUSHROOM 4.8mm/5.8mm (2-4-5-6-8-10)

KITS & SETS

9016-CISK IMPLANT SURGERY KIT 18PCS



1005E EXP23 6EZ



1105E CP15UNC



2037 BONE GAUGE



3506T PERIOSTEALS MOLT9



3509T PERIOSTEALS PR3



4107 LUCAS87



4614 BONE RIDGE MAPPING CALIPER



5301C HEMOSTATIC FORCEPS
BABY MOSQUITO 10cm CVD



2094
BONE MIXING
BOWL



5505B CRILE-WOOD
BLACK EDITION STR 15cm



5911S GOLDMAN-FOX
SUPER CUT STR 12.5cm



6420 UNIVERSITY OF MINNESOTA
RETRACTOR



6501A ADSON CVD TEETH 12cm



6517 LONDON-COLLEGE



6902 3G FLAT SCALPEL HANDLE



6962 #4 MIRROR WITH 6EZ S/E
MIRROR HANDLE (6911 & 6951)



7018N SYRINGE 1.8ml



7512 X CASSETTE
SURGICAL/UTILITY
275x205x35mm

KITS & SETS

9017-CLRA


CROWN LENGTHENING & RIDGE AUGMENTATION KIT 15PCS



1505 YG7/8



1901 PR1/2



2011 OCH2



3509T PERIOSTEALS PR3



3516T PERIOSTEALS BUSER



5303C HEMOSTATIC FORCEPS
HALSTED MOSQUITO 12.5cm CVD




5514B CASTRO STR T/C 18cm
BLACK EDITION



6920 6EZ



6511 SUTURE CORN TISSUE FORCEPS
CVD 16cm



6736 SURGICAL SUCTION
ANDREW PYCHON 23cm



6911 MIRROR HANDLE 6EZ C/S



6951 MIRROR #5 24.8mm
C/S FRONT SURFACE



6525C GERALD MICRO 18cm



7018N SYRINGE 1.8ml



7512 X CASSETTE
SURGICAL/UTILITY
275x205x35mm

KITS & SETS


9018-EBER BERNARD SET




3747 BERNARD 3 LEFT



3748 BERNARD 2 STR



3749 BERNARD CVD



3749IN BERNARD IN



3749OUT BERNARD OUT



3750 BERNARD 4 RIGHT

KITS & SETS

9019-ELUL LUXATOR LARGE HANDLE SET




3751 2mm STR



3752 3mm STR



3753 5mm STR



3754 2mm CVD



3755 3mm CVD



3756 5mm CVD



3785 4mm STR



3786 4mm CVD

KITS & SETS

9020-EXAM EXAM KIT 6PCS



1007E EXP TU17/23 6EZ



1018E ODU 11/12 6EZ




1104E CP12 (3-6-9-12mm)



6517 LONDON-COLLEGE



6962 MIRROR #5 WITH 6EZ HANDLE (6951+6911)



7503S CASSETTE DOUBLE HINGE
6 HOLD SINGLE TALL RACK
90x203x31mm

KITS & SETS


9021-1MCH OSTEOTOME SET 5PCS




2057 OSTEOTOME 4mm STR (7-10-13-15-18mm)



2058 OSTEOTOME 4mm CVD (7-10-13-15-18mm)



2059 OSTEOTOME 7.5mm STR (7-10-13-15-18mm)



2060 OSTEOTOME 7.5mm CVD (7-10-13-15-18mm)



7501 CASSETTE 5 HOLD
60x180x17.5mm

KITS & SETS

9024-IMPL IMPLANT KIT 31 PCS



1105E CP15UNC



1923T PT2 TITANIUM



1924T PT1 TITANIUM



2033N BONE GRAFT PACKER 3.3/4mm



2034S BONE GRAFT PACKER 3/5mm SERRATED



2035T BONE CARRIER TITANIUM



2060 OSTEOTOME 7.5mm CVD (7-10-13-15-18mm)



3506T PERIOSTEALS MOLT9



3751 2mm STR



4103 LUCAS86



4111 MOLT2/4



4618 RULER 15cm 6in



6420 UNIVERSITY OF MINNESOTA RETRACTOR



6428 BISHOP RETRACTOR



6502 ADSON STR TEETH 1x2 12cm



6717 4mm



6905 5A ROUND STR



6911 MIRROR HANDLE 6EZ C/S




6951D MIRROR #5 DOUBLE-SIDED 24.8mm C/S



7018N SYRINGE 1.8ml



4223M MEAD MALLET MINI



5505TN MICRO CRILLE-WOOD STR T/C 15cm

5920 IRIS CVD 11.5cm



2094 BONE MIXING BOWL



2071 KERRISON FORCEP 3mm



7505-1 CASSETTE CLIPS SCISSOR

4625 BONE RIDGE MAPPING CALIPER IMPLANT



7505-2 CASSETTE CLIPS SUCTION



7512 X CASSETTE SURGICAL/UTILITY 275x205x35mm

KITS & SETS

9025-PEEX PERIO EXAM KIT



1006E EXD5 6EZ



1018E ODU 11/12 6EZ



1105E CP15UNC



1125E CPQN (3-6-9-12mm)




6517 LONDON-COLLEGE



6911 6EZ C/S



6951 MIRROR #5 24.8mm C/S FRONT SURFACE



7503S CASSETTE DOUBLE HINGE
6 HOLD SINGLE TALL RACK
90x203x31mm

KITS & SETS

9026-PERA PERIODONTAL SURGERY KIT



1018E ODU 11/12 6EZ



1104E CP12 (3-6-9-12mm)



1125E CPQN (3-6-9-12mm)



1801 K15/16



1812 KOR1/2



1821 F3/4S



2004 C36/37



2010 OCH1



2011 OCH2



3509T PERIOSTEALS PR3



4111 MOLT2/4



6502B ADSON BROWN 12cm



6517 LONDON-COLLEGE



6717 4mm



6911 MIRROR HANDLE 6EZ C/S



6951D MIRROR #5 24.8mm C/S



7502D CASSETTE DOUBLE HINGE
16 HOLD DOUBLE RACK
114x203x31mm

KITS & SETS

9027-PROP LARGE PROPHY KIT 15PCS



1018E ODU 11/12 6EZ



1132E EXP23/CP12



7423

SHARPENING STONE
INDIA SLIP



1201 H6/7



1208 204S



1416 COL13/14



7424

SHARPENING STONE
ARKANSAS FLAT



1601 GR1/2



1604 GR7/8



1606 GR11/12



1607 GR13/14




6517 LONDON-COLLEGE



6911 MIRROR HANDLE 6EZ C/S



6951D MIRROR #5 24.8mm C/S



7502D CASSETTE DOUBLE HINGE
16 HOLD DOUBLE RACK
114x203x31mm

KITS & SETS

9028-SCGS SCALER & CURETTE SET



1201 H6/7



1208 204S



1416 COL13/14



1606 GR11/12



1607 GR13/14



7501 CASSETTE 5 HOLD
60x180x17.5mm

KITS & SETS

9029-SISC IMPLANT SCALER & CURETTE SET



1201T H6/7 TITANIUM



1208T 204S TITANIUM



1416T COL13/14 TITANIUM



1606T GR11/12 TITANIUM



1607T GR13/14 TITANIUM



7501 CASSETTE 5 HOLD
60x180x17.5mm

KITS & SETS

9030-SLEK EXTRACTION KIT



4974
F1 UPPER
UNIVERSAL



4975
F2 LOWER
UNIVERSAL



4978
F5 LOWER
ANTERIOR



4980C
F7C LOWER
MOLAR



4984
F11 UPPER
MOLAR L




4985
F12 UPPER
MOLAR R



4994
46XL LOWER
ROOT



4997
51XL UPPER
ROOT



7511 CASSETTE
SURGICAL/UTILITY
174x203x31mm

KITS & SETS

9032-BBG BONE GRAFT & PACKER SET



2033A BONE GRAFT PACKER 1.0/1.4mm (8-11-13-15-18mm) TITANIUM



2033B BONE GRAFT PACKER 1.7/2.0mm (8-11-13-15-18mm) TITANIUM



2033C BONE GRAFT PACKER 2.0/2.7mm (8-11-13-15-18mm) TITANIUM



2033NT BONE GRAFT PACKER 3.3/4mm TITANIUM



2035T BONE CARRIER TITANIUM



2038T BONE GRAFT/CARRIER 2-4-6-8-10mm TITANIUM

